



**EXCLUSIVE TO STEALTH PRODUCTS**

No matter the season or the situation, the patented and award-winning CoolCore fabric wicks moisture away from your body, helping keep you dry and comfortable. The chemical-free blends of the fabric also help with moisture transportation and moisture evaporation, making our cushions cool better and last longer.

**STEALTH PRODUCTS CUSHION SPECIFICATIONS**

Cushion Name	Weight Capacity	Weight of Cushion	Cushion Height	Width Range	Depth Range
Simplicity G	204 kgs	0.9 kgs	2"	16" - 24"	16" - 20"
Zen SP	204 kgs	1.5 kgs	3"	10" - 24"	10" - 24"
Glacial SP	181 kgs	1.1 kgs	2.5"	10" - 24"	10" - 24"
Premiere P	204 kgs	1.5 kgs	3.5"	16" - 22"	16" - 22"
Spectrum Foam SPP	294 kgs	1.3 kgs	3"	16" - 32"	16" - 28"
Spectrum Gel SPP	226 kgs	1.5 kgs	3"	14" - 24"	16" - 24"
Solution SPP	226 kgs	2.4 kgs	3"	10" - 24"	10" - 22"
Tru Comfort 2 SPP	204 kgs	1.8 kgs	4.5"	16" - 22"	16" - 24"
Essence SPP	226 kgs	1.5 - 2.2 kgs	3.5"	16" - 22"	16" - 24"

**SOLUTIONS FOR LIFE**

Distributed by  
**aidacare™**  
 HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au



### CUSHION COLOUR GUIDE

To make your cushion selection process easier, Stealth Products Cushions are marked with a colour coded patch. The colour of the patch identifies the function.



### SIMPLICITY G

#### Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 16" - 24"W x 16" - 20"D



### ZEN SP

#### Features:

- Slightly Contoured, high-density molded foam base
- FRG Comfort foam topper for skin protection
- Standard Visco topper (exclusive FRG SunMate) for added comfort
- FRG SunMate is an all natural flame retardant material.
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- Optional Urethane liner
- 4-Way Stretch Poly Liner optional
- Sizes available - 10"-24"W x 10"-24"D



### GLACIAL SP

#### Features:

- Contoured, high-density molded foam base
- Soft gel top layer for skin protection
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- 4-Way Stretch Poly or Urethane liner optional
- Sizes available - 10"-24" W x 10"-24"D



### PREMIERE P

#### Features:

- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base
- Built-in positioning platform for off loading
- Standard Visco topper (exclusive FRG SunMate) for added comfort
- FRG SunMate is an all natural flame retardant material.
- Coccyx and sacral relief
- Sizes available - 16"-20"W x 16"-20"D 18" x 22", 20" x 22" 22" x 22" Coming Soon
- CoolCore® technology cover standard.
- Optional Urethane liner
- Optional 4-Way Stretch Poly liner



### SPECTRUM FOAM SPP

#### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 16" - 32"W x 16" - 28"D



### SPECTRUM GEL SPP

#### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 14" - 24"W x 16" - 24"D



### SOLUTION SPP

#### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 22"D
- 4-Way Stretch Poly liner optional



### TRU-COMFORT 2 SPP

#### Features:

- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 16" - 24"D
- 4-Way Stretch Poly liner optional



### ESSENCE SPP

#### Features:

- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base.
- Coccyx and sacral relief.
- CoolCore® technology cover standard.
- 4-Way Stretch Poly liner optional.
- Sizes available - 16"-22"W x 16"-24"D
- Adjustable contour wedge package provides better positioning customization for the pelvic region. Wedge height from 0.5" - 2" with the wedge kit.