

TabTimer™ helps to
keep people on-time™

TabTimer® automatic pill dispensers are widely used throughout Australia, by both social and health care professionals, as a component of a care package for vulnerable patients who have difficulty managing their medicines.

Quality Care Pharmacy Program (QCPP)

Many Pharmacists fill TabTimer® Automatic Pill Dispensers nationwide. The device meets the requirements of the Quality Care Pharmacy Program (QCPP) when used with tamper evident stickers and labelling.

Filler Tray & Tamper Evident Stickers

A pharmacy filler tray and tamper evident stickers are available to assist pharmacy filling.



Careousel Dispenser
Accessories

Browse all filler trays, tamper
evident stickers, and accessories

ANTIMICROBIAL PROTECTION BUILT IN



99.9%
EFFECTIVE

Our Careousel Automatic Pill Dispensers are made with silver antimicrobial technology, laboratory tested to reduce bacteria, mould, and viruses.

1300 822 846

aidacare™
HEALTHCARE EQUIPMENT

Our promise to you

Aidacare's vision is to be Australia's leading healthcare equipment provider, delivering dignity and mobility through quality products, trusted advice and ongoing support.

Aidacare staff are trained and skilled professionals who respond sensitively and expertly to your needs.

Get in touch

1300 822 846



TabTimer Products

Browse all TabTimer
products on our website



Automatic Pill Dispensers



Meets requirements of the
Quality Care Pharmacy Program*

MK3 Careousel

Automatic Pill Dispenser



This device releases the correct dose of medication at the schedule time, designed for safe and reliable medicine administration.

- Up to 24 programmable daily alarms
- 29 compartments for organised storage
- Operates on four included AA batteries

DLG275000

Advance Careousel

Automatic Pill Dispenser



A clear LCD display shows the next-dose time, remaining doses, and battery level to ensure care-friendly monitoring.

- Shutter cover provides extra security
- 24 programmable alarms per day
- 28 compartments for organised storage
- Operates on four included AA batteries

DLG275010

E-Pill Medtime Safe

Tamper Resistant
Automatic Pill Dispenser



Designed for security, this dispenser is suitable for both home and institutional environments managing high-risk medication regimens.

- Robust metal case with multiple independent locks
- 28 programmable alarms per day
- 28 compartments for organised storage
- Operates on four included AA batteries

DLG275020

Accessories



Spare Tray

Pre-fill while main unit is in use

MK3 DLG275100
ADVANCE DLG275110



Tamper Evident Stickers

Ensures medication safety and meets QCPP requirements

DLG275090



Filler Tray

Enables pharmacy pre-filling

MK3 DLG275190
ADVANCE DLG275200



Tipper

Smooth tilt design for low-effort pill access, ideal for limited hand mobility or vision impairment

DLG275080



Frosted Lid

Allows carers to have some visibility of contents

MK3 DLG275160
ADVANCE DLG275170



Pharmacy Label

Meets labelling requirements

Label DLG275140
Holder DLG275150