

WheelAir Product Selection Tool

With a range of WheelAir products available, it is important to choose the right WheelAir package for you. Depending on your wheelchair, it is possible to use both the WA Cushion Cover and Slingback simultaneously. We recommend using the selection tool to decide which product best suits you.

Product Line	Wcare	Wgo	
Product	WA System	WA Cushion Cover	WA Slingback
Description	An integrated modular solution in partnership with our Recognised Partners available in all types of seating systems, such as power chairs and custom seating solutions.	A ready-to-go cover transforming any cushion into a microclimate cushion for optimal temperature and humidity control of the seat surface.	A ready-to-go back support for wheelchairs with velcro straps at the back offering full temperature control for the active user.
User Profile	Struggles with temperature regulation, sweat and/or humidity control	Prone to skin damage at the seat surface	Struggles with temperature regulation
Location	Seat, backrest, sleeping systems	Seat cushion	Backrest support
Channel positioning	Fully customisable	Underneath cushion	At seat surface
Wheelchair fit	All	All, requires breathable cushion	Wheelchairs with velcro straps
Temperature control	++++	+++	++++
Sweat control	++++	++	++++
Humidity control	++++	++++	++++
Airflow sensation	+++	++	++++
Customisation	++++	+	++
Quick installation	++	++++	++++