

## The Perfect Powered Wheelchair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



## Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.\*



When not in use, the Model F can be conveniently stored in a closet or an out-of-the-way corner of the room.



\*Model F is not to be used as a seat while it is being transported. The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.





### Lightweight

Weighing only 24kg, the Model F is one of the lightest power chairs on the market today.\*\*







Top Speed

Range

Obstacle Clearance



780mm turning radius for tight turns

## Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.







Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

# Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



<sup>\*\*</sup>Weight is 26.7kg, if battery is not removed. Battery weighs 2.7kg.

Note: Driving range of 20kms will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

#### Customisation

Choose from five colours and a great selection of accessories to match your personal style and enhance your driving experience.

#### Colour options



Light Blue

WCA870740

WCA870720

Red WCA870730

Controller can be positioned on either the left or right side.





Accessories

Smart Key WCA870700



Cane Holder WCA870710



Lap Belt WCA870420



Battery\* WCA870400



Battery Charger\* WCA870410



**Bottom Basket** MSP492510

\*Included as standard

## Specifications

Product Code	PWS655500
Driving Range	Up to 20km
Foldable	Yes
Turning Radius	780mm
Seat Sizes Available	450mm / 18"
Device Width	555-605mm
Device Weight	26.7kg
App Enabled	Yes, with remote control
Max Weight Capacity	Up to 113kg
Battery	Lithium-Ion (25.3V 10Ah)

Maximum Speed***	8km/h
Max Incline	10°
Obstacle Clearance	35mm
Device Length	935mm (465mm when folded)
Charging Time	~ 5 hours
Air Travel Approved	Yes
Swappable Arm Covers	5 Colours
Water Resistance	IPX4
Width Between Armrests	Adjustable: 2 positions (400-450mm)

<sup>\*\*\*</sup>Acceleration and/or decleration may change when driving on steep slopes.

<sup>2.</sup> The device is designed for use on ordinary paved roads or indoor use.











<sup>1.</sup> This device is intended to provide mobility to individuals who are comfortable with operating a powered wheelchair. Persons who find it difficult to use a powered wheelchair should refrain from using the device.