


4H400E Shower Trolley - Electric - 730mmW x 1890mmL - 150kg

Product Code	Product Description
BTS147100	4H400E Shower Trolley - Electric - 730mmW x 1890mmL - 150kg ^



Shower Trolley w/TiMotion Electronics (Current version)

No.	Item Number	Item Description
1	BTS147110	Liner - Shower trolley - Standard - 4Healthcare & Kerry
2	GSP005750	Side Rail - 4H400E - 4HC Shower Trolley
3	GSP005752	End Rail - 4H400E - 4HC Shower Trolley
4	GSP005754	Round Buffer - (Wall Bumper) - 100mm - 4H400E - 4HC Shower Trolley
5	GSP005722	Castor - Braking - Tente - 125mm - suits 4H400E - 4HC Shower Trolley
6	GSP005723	Castor - Directional - Tente - 125mm - suits 4H400E - 4HC Shower Trolley
7	GSP005756	Shower Pillow - 4H400E - 4HC Shower Trolley
8	GSP005758	Lift Coulmn - TiMotion - 4H400E - 4HC Shower Trolley
9	MBT401050	Battery - TiMotion
9A	GSP005764	Battery box - 4H400E - 4HC Shower Trolley
10	GSP005760	Handset - 4H400E - 4HC Shower Trolley
11	GSP005762	Control Box - 4H400E - 4HC Shower Trolley
12	BTS147130	Hose Kit - 4HC Shower Trolley

Shower Trolley w/Linak Electronics (Older version)

No.	Item Number	Item Description
8	GSP874332	Lifting Column
9	GSP342744	Battery - NEW STYLE - 4H400E Shower Trolley
10	GSP342735	Handset
11	GSP342742	Control Box

This manual is for internal use only.