

## Aspire Nova Standing Lifter Spare Parts Manual

---




Item Number	Item Description
LSS390652	Aspire Nova - Standing Lifter - 250kg


**This manual is for internal use only**



Spare Parts Diagram




# 1 Base and Castors



Ref	Item Number	Item Description	Image	Type
				




Ref	Item Number	Item Description	Image	Type
1	GSP242800	Front Castor - w/Mounting Hardware - Aspire Nova Standing / Aura Flor Lifter		Spare Part
2	GSP242802	Rear Castor - w/Mounting Hardware - Aspire Nova Standing / Aura Floor Lifter		Spare Part
3	GSP242804	Base Frame - Complete - w/Castors - w/o Electronics - Aspire Nova Standing Lifter		Spare Part

1 Base and Castors				
Ref	Item Number	Item Description	Image	Type
4	GSP242806	Plastic Covers - Legs - Pair - Aspire Nova Standing Lifter		Spare Part
5	GSP242808	Rubber Gromet - on Base Frame - Aspire Nova Standing / Aura Flor Lifter		Spare Part
6	GSP242810	Foot Board Platform - Aspire Nova Standing Lifter		Spare Part
6	GSP242812	Rubber Liner - Foot Board - Aspire Nova Standing Lifter		Spare Part
7	GSP242814	Hardware - Legs Pivot Point - Aspire Nova Standing / Aura Flor Lifter		Spare Part
8	GSP242816	Plastic Gromet - Base to Mast - Aspire Nova Standing Lifter		Spare Part


## 2 Mast, Boom, Handle and Knee Pad Assembly

Ref	Item Number	Item Description	Image	Type
				

Ref	Item Number	Item Description	Image	Type
15	GSP242824	Mast Frame w/Pivot Insert (on Top End) - Aspire Nova Standing Lifter		Spare Part
16	GSP242826	Mounting Hardware only - Hinge Assembly - Aspire Nova Standing / Aura Floor Lifter		Spare Part
17	GSP242828	Boom Frame / Sling Bar - w/Hook Insert - (Incl. Pivot Insert) - Aspire Nova Standing Lifter		Spare Part

2 Mast, Boom, Handle and Knee Pad Assembly				
Ref	Item Number	Item Description	Image	Type
17	LSS390660	Boom Frame / Sling Bar - w/Clip Insert - (Incl. Pivot Insert) - Aspire Nova Standing Lifter		Accessory
18	GSP242832	End Cap - Boom Frame / Sling Bar - w/Clip Insert - Aspire Nova Standing Lifter		Spare Part
19	GSP242834	Handle Grip - Pair - Aspire Nova Standing Lifter		Spare Part
20	GSP242836	Handle Bar - w/Mounting Hardware - Aspire Nova Standing Lifter		Spare Part
21	GSP242838	Knee Pad Assembly - Complete - w/Bracket Frame, Pad - Aspire Nova Standing Lifter		Spare Part
22	GSP242840	Neoprene Calf Strap - Aspire Nova Standing Lifter		Spare Part






## 2 Mast, Boom, Handle and Knee Pad Assembly

Ref	Item Number	Item Description	Image	Type
24	GSP242844	Hand Wheel - for Knee Pad Adjustment - Aspire Nova Standing Lifter		Spare Part



### 3 Electronics and Hardware

Ref	Item Number	Item Description	Image	Type

Ref	Item Number	Item Description	Image	Type
32	GSP242852	Lifting Actuator - Aspire Nova Standing Lifter		Spare Part
32A	GSP242854	Mounting Hardware - Top and Bottom - w/Cap - Lifting Actuator - Aspire Nova Standing / Aura Floor Lifter		Spare Part
32B	GSP242856	Plastic Cap - (Bolt Cover) - Aspire Nova Standing / Aura Floor Lifter		Spare Part

3 Electronics and Hardware				
Ref	Item Number	Item Description	Image	Type
33	GSP242858	Leg Actuator (Only) - Aspire Nova Standing / Aura Floor Lifter		Spare Part
33A	GSP242860	Mounting Hardware - w/Bushing - Leg Actuator - Aspire Nova Standing / Aura Floor Lifter		Spare Part
34	GSP242862	Control Box - Aspire Nova Standing / Aura Floor Lifter		Spare Part
36	MBT401050	Battery -Aspire Nova Standing / Aura Floor Lifter		Spare Part
37	MBT401055	Charging Cradle - w/Transformer Cable - Aspire Nova Standing / Aura Floor Lifter		Spare Part
37A	MBT401060	Power Transformer - w/Cable - for Charging Cradle - Aspire Nova Standing / Aura Floor Lifter		Spare Part

**3 Electronics and Hardware**

Ref	Item Number	Item Description	Image	Type
38	GSP242872	Handset - 4 Button - Aspire Nova Standing / Aura Floor Lifter		Spare Part
39	GSP242874	Splitter Cable - 3 Channel - (Wiring Harness) - Aspire Nova Standing / Aura Floor Lifter		Spare Part