


Aidacare Aspire A320 Bariatric Lifter (LSS390630)

Spare Parts Manual











Aidacare Aspire A320 Bariatric Lifter – Exploded Parts Diagram








Aidacare Product Code	Aidacare Product Description
LSS390630	Patient Lifter - Aidacare Aspire A320 - Alum - Electric Leg - 320kg

Spare Parts List



Ref #	Item Number	Item Description	Type
1, 2	VARIOUS	Castor Assembly - Tente	Spare Part
<p>PLEASE REFER TO FILE NAMED 'TENITE CASTOR ASSEMBLY' ON INTRANET FOR CORRECT PART NUMBERS AND CONSUMABLES TO USE FOR CASTORS</p>			
3	GSP390933	Metal Castor Bracket – Rear & Front	Spare Part
			
4	GSP390930	Leg Hinge Bolt Assembly	Spare Part
			
5	GSP390911	Plastic Cap - Bolt Cover	Spare Part
			
6	GSP390912	Plastic Shroud for Base of Mast	Spare Part
			





7	GSP390915	Plastic Castor Cover – Front – PAIR	Spare Part
			
8	GSP390918	Plastic Castor Cover – Rear – PAIR	Spare Part
			
9	GSP390920	Plastic Cover – Leg Gate Actuator	Spare Part
			
9a	GSP390677	Electric Leg kit Bracket	Spare Part
			

10	GSP390908	Sling / Spreader Bar	Spare Part
10d	GSP390868	Handling Velcro Pads – for Sling / Spreader Bar – Set of 4	Spare Part
			
10a	GSP390976	Sling Bar Retainer Clip	Spare Part
			
10b	GSP390991	Adaptor Pin	Spare Part
			
10c	GSP390993	Link Pin Spreader Bar to Mast Assembly	Spare Part
			
11	GSP390940	Handgrip Padding – PAIR	Accessory
			

12	GSP390691	4 Button Handset	Spare Part
12A	GSP390692	Programming Handset - for Aspire Lifters	
			
13	MBT401200	Battery with Port	Spare Part
			
13	MBT401190	Battery without Port	Spare Part
13a	MBT401220	Battery Charger - Plug Pack	Spare Part
13b	MBT401210	Battery Charger - Wall Mount	Spare Part
13c	MBT401012	Battery Charger - Linak - For Lifters - Including Mounting Brackets & Cable	Spare Part
14	GSP390688	Control Box	Spare Part
			
14a	MSP500906	Bracket – Short – for Control Box & Battery	Spare Part
15a	GSP394097	Cable Straight – for Leg Actuators	Spare Part

			
15b	GSP394172	Cable Curly – for Lifter Actuators	Spare Part
			
15	GSP390687	Actuator - Lifter	Spare Part
			
16	GSP390921	Actuator Bolt Assembly – Upper	Spare Part
17	GSP390941	Actuator Bolt Assembly – Lower	Spare Part
			
18	GSP390679	Actuator – Leg Gate	Spare Part

			
19	GSP390946	Adjustment Lever – Mast	Spare Part
			
20	GSP390953	Replacement Leg - Each	Spare Part
			
21a	GSP390961	Plastic Hinge Covers	Spare Part
			
21b	GSP390669	Plastic Hinge Cover w Slot Cut Out	Spare Part
			
21c	GSP390987	Hinge Assembly Mounting Hardware	Spare Part
			
22	GSP390965	Tool Assembly Kit	Spare Part

				
23	GSP390876	Boom Frame		Spare Part
				
24	GSP390877	Mast Frame		Spare Part
				
25	GSP390996	Handle Bar		Spare Part
				
26	GSP390997	End Caps - Handle Bar		Spare Part
				
30	GSP390790	Replacement Sticker Pack		Spare Part

aidacare
aspire

WARNING!

ALWAYS read the instruction manual if you are unfamiliar with how to use this lifter.

ALWAYS use Aidacare approved slings with this lifter. Incorrect slings may cause serious injury and/or death. For more information, see the warranty. If in doubt, contact Aidacare 1300 133 120.

ALWAYS inspect lifter and sling prior to use.

DO NOT use lifter on uneven surfaces or on surfaces where there is a risk of the lifter could tip over.

DO NOT USE if damage is suspected to or visible on the lifter or sling.

IF IN DOUBT - DO NOT USE THIS LIFTER.

SWL 320KG

WARNING!

aidacare
aspire **A320**

MAINTENANCE AND IMPROVE

Aidacare Pty Ltd
ABN: 60 134 398 888
84 King St,
Moorebank, NSW
Moorebank NSW 2170

CONTACT DETAILS

T: 1300 133 120
F: 62 9618 6111
E: product@aidacare.com.au

SALES

It is recommended that you have this product serviced every 6 months.

aidacare (LIFE) CARE

SMALL 235kg	MID 295kg	XL 395kg
LOCK 395kg	KEVIN LAUREN 395kg	

aidacare
HEALTHCARE EQUIPMENT

1300 133 120 www.aidacare.com.au

SWL 320KG

CAUTION!

ONLY USE ASPIRE Lifter with its own base. Be careful when assembling & keep fingers clear of the mast receiver bracket when inserting the mast.

DO NOT STAND ON COVER

SAFE WORKING LOAD 320KG

DANGER

PINCH POINT