


Active Air 4 Mattress Overlay - Spare Parts Diagram



Product Code	Product Description
PTM639100	Mattress Overlay - Aspire ACTIVE AIR 4 - Single - 19 Cells

No	Product Code	Product Description
1	GSP493040	Top Cover - Aspire Active Air 4 Mattress
2	GSP493090	Base Cover - Active Air 4 Mattress
2-1	GSP493420	Cell Retainer - Active Air 4
3	GSP493570	Mattress Pump Hose & Connector - Aspire Active Air 4 Mattress
4	GSP493830	Pump Hose Sleeve Cover - Aspire Active Air Mattresses
5	GSP493600	CPR Twist Valve
6	GSP493620	CPR Valve Sleeve Cover - Aspire Active Air Mattresses
7	GSP493370	Mattress Cell Manifold - Aspire Active Air 4 Mattress
8	GSP493150	Replacement Cell - Dual Left (1) /Right (1) Connectors - Head/Foot - Aspire Active Air 4 Mattress
9	GSP493140	Replacement Cell - Single Left (1) Connector - Body - Aspire Active Air 4 Mattress
10	GSP493160	Replacement Cell - Single Right (2) Connectors - Body - Aspire Active Air 4 Mattress
11	GSP493150	Replacement Cell - Dual Left (1) /Right (1) Connectors - Head/Foot - Aspire Active Air 4 Mattress

Product Code	Product Description
PTM639010	Carry Bag - Aspire Active Air 4/5/5KS Mattresses

Prior SRN: 1209436885 Button Stud Colour - Grey or Fluorocent Green
 From SRN: 1209436885 Button Colour - Black
 Button Cap Remains Green