

### Spare Parts Breakdown - Aspire Aluminium Commode - Att and Self Propelled



	BTC068010	Shower Commode - Aspire Aluminium - Att Prop - 440mm (seat incl)
	BTC068020	Shower Commode - Aspire Aluminium - Self Prop - 440mm (seat incl)

No.	Product Code	Product Description
1	GSP491210	Anti-tippers
2	GSP491220	Backrest & Mounting Hardware
3	GSP491230	Commode Bowl & Lid
4	GSP491242	Fastening Hand Wheels - Height / Depth / Backrest - ROUND / STAR
5	GSP491250	Front Castor (non-braking) & Mounting Hardware
6	GSP491260	Padded Armrest & Mounting Hardware - Left
6	GSP491270	Padded Armrest & Mounting Hardware - Right
7	GSP491280	Padded Seat & Mounting Hardware
7	GSP491290	Padded Seat Plug
8	GSP491300	Quick Release Axle
9	GSP491310	Rear Brake & Mounting Hardware - Left
9	GSP491320	Rear Brake & Mounting Hardware - Right
10	GSP491330	Rear Castor (braking) & Mounting Hardware
10	GSP491340	Rear Self Propelling Wheel
11	GSP491350	Swing Away Legrest - Left
11	GSP491360	Swing Away Legrest - Right
12	GSP491345	Lock Brake for Swing Away Legrest
13	GSP491346	Footplate Strap
14	GSP491347	Backrest Plug
15	GSP491348	Armrest End Cap