

Aspire Lift Recline Chair DaVinci - Spare Parts Diagram



Product Code	Product Description
CHP227500	Chair Powerlift Recline - Aspire Da Vinci Quattro - Latte - 4 Motor - 159kg
CHP227510	Chair Powerlift Recline - Aspire Da Vinci Quattro - Gunmetal - 4 Motor - 159kg
CHP227520	Chair Powerlift Recline - Aspire Da Vinci Quattro - Fog - 4 Motor - 159kg
CHP227530	Chair Powerlift Recline - Aspire Da Vinci Quattro - Dove - 4 Motor - 159kg

No	Item Number	Product Description
1	GSP510385	Hand Control - Aspire Lift Recline Chair DaVinci
2	GSP510425	Transformer - Aspire Lift Recline Chair DaVinci
3	GSP510085	Actuator - Hi/Lo - Aspire Lift Recline Chair DaVinci
4	GSP510090	Actuator - Backrest - Aspire Lift Recline Chair DaVinci
5	GSP510550	Frame - Base Chair - Aspire Lift Recline Chair DaVinci
6	GSP510560	Frame - Headrest and Lumbar - Aspire Lift Recline Chair DaVinci
7	GSP510125	Cable - Aspire Lift Recline Chair DaVinci - for Transformer to Control Box
8	GSP510426	Power Cord - Aspire Lift Recline Chair DaVinci
9	GSP510045	Actuator - Headrest - Aspire Lift Recline Chair DaVinci
10	GSP510055	Actuator - Lumbar - Aspire Lift Recline Chair DaVinci
11	GSP510237	Control Box - Aspire Lift Recline Chair DaVinci
12	GSP510570	Bracket - Aspire Lift Recline Chair DaVinci - for Pivot attachment on seat and Backrest
13	GSP510333	Cover - Headrest/Armrest - Aspire Lift Recline Chair DaVinci - Latte
13	GSP510334	Cover - Headrest/Armrest - Aspire Lift Recline Chair DaVinci - Gunmetal
13	GSP510335	Cover - Headrest/Armrest - Aspire Lift Recline Chair DaVinci - Fog
13	GSP510336	Cover - Headrest/Armrest - Aspire Lift Recline Chair DaVinci - Dove
14	GSP511845	Side Wings - Aspire DaVinci Lift Recline Chair - Left - LATTE
15	GSP511850	Side Wings - Aspire DaVinci Lift Recline Chair - Right - LATTE
16	GSP511855	Side Wings - Aspire DaVinci Lift Recline Chair - Left - GUNMETAL
17	GSP511860	Side Wings - Aspire DaVinci Lift Recline Chair - Right - GUNMETAL