


Spare Parts Diagram

BEB046030


Aspire LIFESTYLE Community Bed w/Transport Bracket – Gunmetal Grey – Single

BEB046040

Aspire LIFESTYLE Community Bed w/Transport Bracket – Gunmetal Grey – King Single



47



More Spare Parts



2.1



42



45



43



44



46



11.1




47




3.1



5.1



HC_PO number_000001 (sequential until end of range)




40



5.2




HC_PO number_000001 (sequential until end of range)



41

Accessory



7



1

Item Code	Item Description
BEB046030	Bed - Aspire LIFESTYLE Community Bed w/Transport Bracket - Gunmetal Grey - Single - 220kg
BEB046040	Bed - Aspire LIFESTYLE Community Bed w/Transport Bracket - Gunmetal Grey - King Single - 220kg

Bed Frame (Page 1)				
No.	QTY	Item Code	Item Description	Type
1	2	GSP392312	Wood Panels only w/screws - for Bed Ends - Aspire LIFESTYLE Community Care Bed – Single	Spare Parts
1	2	GSP392314	Wood Panels only w/screws - for Bed Ends - Aspire LIFESTYLE Community Care Bed - King Single	Spare Parts
3	4	GSP392310	Mattress Retainers - Side Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
3.1	1	GSP392315	Mattress Retainer - Foot End - Aspire LIFESTYLE Community Care Bed	Spare Parts
4	2	GSP392240	Frame Connection Bar - for Headend / Footend Platform - Aspire LIFESTYLE Community Care Bed	Spare Parts
5	2	GSP392210	Recline Linkages - For Backrest Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
6	1	GSP392250	Platform Assembly - Knee Bend - Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
6	1	GSP392255	Platform Assembly - Knee Bend - King Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
7	1	GSP392200	Platform Assembly - Backrest - Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
7	1	GSP392205	Platform Assembly - Backrest - King Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
8	1			Spare Parts
9	1	GSP392260	Platform Assembly - Foot End - Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
9	1	GSP392265	Platform Assembly - Foot End - King Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
10	1	GSP392290	Frame Connection Bar - for Mattress Extension Platform - Aspire LIFESTYLE Community Care Bed	Spare Parts
11	1	GSP392290	Frame Connection Bar - for Mattress Extension Platform - Aspire LIFESTYLE Community Care Bed	Spare Parts
11.1	1	GSP392300	Mattress Extension Platform - King Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
11.1	1	GSP392295	Mattress Extension Platform - Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
14	2	GSP392210	Recline Linkages - For Backrest Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
15	2	GSP392215	Rubber Retention Strap and Pin - Aspire LIFESTYLE Community Care Bed	Spare Parts
17	1	GSP392270	Rachet Assembly Kit - For Knee Bend - Aspire LIFESTYLE Community Care Bed	Spare Parts
18	1	GSP392270	Rachet Assembly Kit - For Knee Bend - Aspire LIFESTYLE Community Care Bed	Spare Parts
19	4	GSP392215	Rubber Retention Strap and Pin - Aspire LIFESTYLE Community Care Bed	Spare Parts
20	3	GSP392026	Actuator Mounting Hardware (Quick Release Clip) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts
21	3	GSP392026	Actuator Mounting Hardware (Quick Release Clip) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts
22	2	GSP392290	Frame Connection Bar - for Mattress Extension Platform - Aspire LIFESTYLE Community Care Bed	Spare Parts
23	4	GSP392215	Rubber Retention Strap and Pin - Aspire LIFESTYLE Community Care Bed	Spare Parts
24	2	GSP392230	Sliding Block - for Backrest - Aspire MAXI LIFESTYLE / LIFESTYLE Community Care Bed	Spare Parts
25	4	GSP392247	Cable Retainer Clips - Aspire MAXI LIFESTYLE Community Care Bed	Spare Parts
26	4	GSP392275	Rubber Caps - White - Aspire MAXI LIFESTYLE / LIFESTYLE Community Care Bed	Spare Parts

27	2	GSP392225	Plastic Bush Insert - for Self Help / IV Pole - Aspire LIFESTYLE Community Care Bed	Spare Parts
28	4	GSP392220	End Cap Kit - For Mattress Platform - Aspire LIFESTYLE Community Bed	Spare Parts

29	7	GSP392222	End Cap Kit - For Mattress Frame - Aspire LIFESTYLE Community Bed	Spare Parts
31	4	GSP392220	End Cap Kit - For Mattress Platform - Aspire LIFESTYLE Community Bed	Spare Parts
32	2	#N/A	Common Hardware	Spare Parts
33	6	GSP392235	Hand Wheel Bolt - for Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
34	2	GSP392215	Rubber Retention Strap and Pin - Aspire LIFESTYLE Community Care Bed	Spare Parts
35	6	GSP392210	Recline Linkages - For Backrest Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
36	4	#N/A	Common Hardware	Spare Parts
37	2	GSP392210	Recline Linkages - For Backrest Frame - Aspire LIFESTYLE Community Care Bed	Spare Parts
38	6	GSP392215	Rubber Retention Strap and Pin - Aspire LIFESTYLE Community Care Bed	Spare Parts
39	14	#N/A	Common Hardware	Spare Parts
40	1	GSP392305	Sticker Pack - Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
41	1	GSP392306	Sticker Pack - King Single - Aspire LIFESTYLE Community Care Bed	Spare Parts
46	1	GSP392245	Cable Retainer Support Bush - Aspire LIFESTYLE Community Care Bed	Spare Parts
47	1	GSP392280	Locking Pins - for Mattress Extension Platform - Aspire LIFESTYLE Community Care Bed	Spare Parts
1	1	GSP392130	Foot Actuator - Linak - Aspire LIFESTYLE / COMMUNITY CARE Bed	Electronics
2	1	GSP392125	Back Actuator - Linak - Aspire LIFESTYLE / COMMUNITY CARE Bed	Electronics
2.1	1	GSP392135	Cable - Back / Foot Actuator - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
42	1	GSP392145	Control Box - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
43	1	GSP392155	Handset - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
44	1	GSP392150	Mains Power Cable - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
45	1	GSP392160	Plastic Key - for Handset - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
47	1	GSP392027	Bed End Complete - w/Frame, Boards, Castors and Actuator - Aspire LIFESTYLE Community Bed - Single	Spare Parts
47	1	GSP392028	Bed End Complete - w/Frame, Boards, Castors and Actuator - Aspire LIFESTYLE Community Bed – King Single	Spare Parts

Head / Foot Boards (Page 2)				
No.	QTY	Item Code	Item Description	Type
1	2	BEA029850	Bed End Shrouds - Laminate - Single - Aspire LIFESTYLE Community Bed	Spare Parts
1	2	BEA029855	Bed End Shrouds - Laminate - King Single - Aspire LIFESTYLE Community Bed	Spare Parts
2	2			Spare Parts
3	2			Spare Parts
4	2	GSP392025	Actuator Mounting Hardware (Quick Release Pin) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts

5	2	GSP392120	Head / Foot End Actuator - Hi Lo - Linak - Aspire LIFESTYLE / COMMUNITY CARE Bed	Electronics
5.1	1	GSP392140	Cable - Foot End Actuator - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
5.2	1	GSP392142	Cable - Head End Actuator - Linak - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
7	2	GSP392025	Actuator Mounting Hardware (Quick Release Pin) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts
8	4	#N/A	Common Hardware	Spare Parts

9	8	#N/A	Common Hardware	Spare Parts
10	4	GSP392335	Spring Locking Pin - for Side Rails - Aspire LIFESTYLE Community Bed	Spare Parts
11	2	GSP392026	Actuator Mounting Hardware (Quick Release Clip) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts
12	2	GSP392026	Actuator Mounting Hardware (Quick Release Clip) - Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Spare Parts
13	4	GSP392335	Spring Locking Pin - for Side Rails - Aspire LIFESTYLE Community Bed	Spare Parts
14	4			Spare Parts
15	4	GSP392320	Locking Safety Circlip - for Folding - Aspire LIFESTYLE Community Bed	Spare Parts
16	8	#N/A	Not Available	Spare Parts
17	4	GSP392335	Spring Locking Pin - for Side Rails - Aspire LIFESTYLE Community Bed	Spare Parts
18	4	GSP392219	End Cap Kit - For Head / Foot Boards - Aspire LIFESTYLE Community Bed	Spare Parts
19	4	#N/A	Common Hardware	Spare Parts
20	8	#N/A	Common Hardware	Spare Parts
21	4	GSP392325	Castor - Total Locking - Aspire LIFESTYLE Community Bed	Spare Parts
22	4	GSP392330	Handwheel Bolt - for Side Rails - Aspire LIFESTYLE Community Bed	Spare Parts
23	4	GSP392325	Castor - Total Locking - Aspire LIFESTYLE Community Bed	Spare Parts
24	4	#N/A	Common Hardware	Spare Parts
25	4	GSP392325	Castor - Total Locking - Aspire LIFESTYLE Community Bed	Spare Parts

Accessory (Page 4 & Aidacare Brochure's)				
No.	QTY	Item Code	Item Description	Type
1	1	GSP392218	End Cap Kit - For Transport Bracket - Aspire LIFESTYLE Community Bed	Spare Parts
2	1	BEA016596	Bed Self Help Pole - Aspire LIFESTYLE Community Care Bed	Accessory
3	2	BEA029850	Bed End Shrouds - Laminate - Single - Aspire LIFESTYLE Community Bed	Accessory
4	2	BEA029855	Bed End Shrouds - Laminate - King Single - Aspire LIFESTYLE Community Bed	Accessory
5	2	BEA029863	Bed Side Rails - Full Length - Laminate - Aspire LIFESTYLE Community Care Bed	Accessory
5A		GSP392045	Full Length Side Rail Mounting Hardware - Full Length Beech Side Rails - Aspire Community/Homecare	Spare Part
6	2	BEA029864	Bed Side Rail Extension Kit - Laminate - Aspire LIFESTYLE Community Bed	Accessory
7	2	BEA029871	Bed Transport Kit - Aspire LIFESTYLE Community Bed	Accessory