


Spare Parts Breakdown - Aspire Lite Transit Wheelchair



MWS442301 Aspire Lite Transit Wheelchair	
1	GSP491000 Backrest Upholstery
2	GSP491005 Seatrest Upholstery
3	GSP491010 Armrest Assembly & Mounting Hardware - LEFT
3	GSP491015 Armrest Assembly & Mounting Hardware - RIGHT
4	GSP491020 Front Wheel Assembly, Fork & Mounting Hardware
5	GSP491025 Seat Belt
6	GSP491070 Rear Wheel
7	GSP491075 Rear Axle & Mounting Hardware - LEFT
7	GSP491080 Rear Axle & Mounting Hardware - RIGHT
8	GSP491085 Rear Brake Assembly & Mounting Hardware - LEFT
8	GSP491090 Rear Brake Assembly & Mounting Hardware - RIGHT
9	GSP491095 Backrest Cane, Brake Handle & Cable - LEFT - GREEN
9	GSP491100 Backrest Cane, Brake Handle & Cable - RIGHT - GREEN
9	GSP491096 Backrest Cane, Brake Handle & Cable - LEFT - CHAMPAGNE
9	GSP491099 Backrest Cane, Brake Handle & Cable - RIGHT - CHAMPAGNE
9a	GSP491102 Mounting Hardware Bracket - Left - for Backrest Cane
9b	GSP491103 Mounting Hardware Bracket - Right - for Backrest Cane
10	GSP491105 Swing Away Footrest LEFT - GREEN <i>(Discontinued, use the existing stock first)</i>
10	GSP491122 Swing Away Footrest LEFT - GREEN <i>(Discontinued, use the existing stock first)</i>
10	GSP491110 Swing Away Footrest RIGHT - GREEN <i>(Discontinued, use the existing stock first)</i>
10	GSP491120 Swing Away Footrest RIGHT - GREEN <i>(Discontinued, use the existing stock first)</i>
10	GSP491106 Swing Away Footrest LEFT - CHAMPAGNE <i>(Discontinued, use the existing stock first)</i>
10	GSP491126 Swing Away Footrest LEFT - CHAMPAGNE <i>(Discontinued, use the existing stock first)</i>
10	GSP491111 Swing Away Footrest RIGHT - CHAMPAGNE <i>(Discontinued, use the existing stock first)</i>
10	GSP491124 Swing Away Footrest RIGHT - CHAMPAGNE <i>(Discontinued, use the existing stock first)</i>
11	GSP491098 Brake cable set