


Aspire Lift Recline Chair MAUI - Spare Parts Diagram



Product Code	Product Description
CHP227110	Chair Powerlift Recline - Aspire Maui - Mocha - Dual Motor - 130kg

No.	Item Number	Item Description
1	GSP510000	Actuator - Backrest - Aspire Lift Recline Chair MAUI
2	GSP510020	Actuator - Footrest - Aspire Lift Recline Chair MAUI
3	GSP510130	Castor - Locking - Aspire Lift Recline Chair MAUI/RAPHAEL/VERONA/RENOIR
4	GSP510180	Castor - Standard - Aspire Lift Recline Chair MAUI/RAPHAEL/VERONA/RENOIR
5	GSP510452	Cover Set - Armrest / Headrest - Aspire Lift Recline Chair MAUI
6	GSP510452	Cover Set - Armrest / Headrest - Aspire Lift Recline Chair MAUI
7	GSP510350	Plastic Foot - Aspire Lift Recline Chair MAUI/MONTANA/RAPHAEL / RENOIR/VERONA
8	GSP510355	Plastic End Cap - Side Frame - Aspire Lift Recline Chair MAUI/MONTANA/RAPHAEL / RENOIR/VERONA
9	GSP510381	Handset - Aspire Lift Recline Chair MAUI - NEW with Aspire Print
10	GSP510404	Mains Power Plug (AU) - Aspire Lift Recline Chair MAUI/MONTANA/VERONA
11	GSP510411	Power Pack / Transformer - Aspire Lift Recline Chair MONTANA/VERONA/MAUI
12	GSP400648	Transformer ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds
13	GSP400649	Battery ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds
14	GSP510430	Seat/Leg Upholstery - Aspire Lift Recline Chair MAUI
15	GSP511240	Base Frame - Aspire Lift Recline Chair MAUI
15.1	GSP511348	Backrest Frame - Aspire Lift Recline Chair MAUI