


No	Code	Description
1	GSP046300	Padded Backrest Assembly - Aspire MAXI Adjustable Bedside Commode - 550mm - Champagne
2	GSP046302	Padded Seat - Aspire MAXI Adjustable Bedside Commode - 550mm - Champagne
3	BTC066755	Shower Commode Seat - Aspire MAXI Adjustable Bedside Commode - Open - 550mm
4	GSP046304	Padded Armrests - Aspire MAXI Adjustable Bedside Commode - PAIR - Champagne
5	GSP145020	Frame only - Aspire MAXI Adjustable OTA & Bedside Commode - 550 to 850mm - 310kg
7	GSP145010	Replacement Leg - Aspire MAXI Shower Chair / Stool / OTA / Commode - Set/4
8	GSP144990	Rubber Ferrule / Tip - Aspire MAXI Shower Chair / Stool / OTA / Commode
9	BTT145105	Plastic Splash Guard only - Aspire - for MAXI Range (BTT145020/BEC045125)
10	BTC081610	Bowl & Lid - Aspire - Classic/Deluxe Bedside/Shower Commodes ^
11	GSP144940	Backrest Adjustment Pin - Aspire MAXI Shower Chair/Bedside Commode
12	GSP144960	End Rubber Cap - Aspire MAXI Shower Chair / Stool / OTA / Commode
13	GSP144970	C-Clip Pin - Aspire MAXI Shower Chair / Stool / OTA / Commode
14	GSP144950	Fastening Hand Wheel - for Seat Frame - Aspire MAXI Shower Chair / Stool / Commode



	NAME	SIGNATURE	DATE	 <small>AIDACARE PTY LTD BUILDING 3A 1 MOOREBANK AVE MOOREBANK NSW 2170 AUSTRALIA</small>
DRAWN	MZ		2024-01-03	
CHECKED				
APPROVED	MATTI		2024-12-16	TITLE SPARE PARTS DIAGRAM BEC045125 - CHAMPAGNE MAXI BEDSIDE COMMODOE
	MATERIAL	FINISH	SIZE	DWG NO.
	N/A	AS MANUFACTURED	A3	D04297
			SCALE	WEIGHT
			1:14	
			SHEET	REV.
			1 of 1	J