


Aspire MAXI Range - Spare Parts Diagram



Product Code	Product Description
BTS118020	Aspire Maxi Shower Chair
BTS119030	Aspire Maxi Shower Stool
BTT145020	Aspire Maxi Over Toilet Aid

No	Code	Spare Part Description	Shower Chair	Shower Stool	Over Toilet Aid
1	GSP145040	Padded Backrest Assembly - w/brackets - 550mm	Yes	No	No
1	GSP145050	Padded Backrest Only - w/o brackets - 550mm	Yes	No	No
2	BTC066750	Padded Toilet Seat - Open Seat - 530mm	No	No	Yes
2	GSP145030	Padded Seat - 550mm	Yes	Yes	No
2(i)	GSP145035	Seat Assembly	Yes	Yes	Yes
3	GSP145060	Padded Armrests - White Vinyl	Yes	Yes	No
3(i)	GSP145020	Frame Only - 550 to 850mm	Yes	Yes	Yes
4	GSP144980	Plastic Bushing for Leg Frame	Yes	Yes	Yes
5	GSP145000	Replacement Leg - Single 6 Holes	Yes	Yes	Yes
5	GSP145010	Replacement Leg - SET / 4	Yes	Yes	Yes
6	GSP144990	Rubber Ferrule/ Tip	Yes	Yes	Yes
7	BTT145105	Splash Guard	No	No	Yes
8	BTC081610	Bowl & Lid	No	No	Yes
9	GSP144940	Backrest Adjustment Pin	Yes	No	No
10	GSP144950	Fastening Handwheel for Seat Frame	Yes	Yes	No
11	GSP144960	End Cap	Yes	Yes	Yes
12	GSP144970	C-Clip Pin	Yes	Yes	Yes