












## **Spare Parts Diagram**

**BEB046020**

**Aspire MAXI Lifestyle Community Bed – Electric Folding – Gunmetal Grey – 385kg**

No.  
QTY





20



40



41




21/22



24a



24a



HC\_PO number \_000001 (sequential until end of range)



64

Product Code	Product Description
BEB046020	Bed Aspire MAXI LIFESTYLE Community Bed Electric Folding Gunmetal Grey 385kg

Item	QTY	Item Code	Item Description	Type
1	4	GSP392310	Mattress Retainers Side Frame Aspire MAXI LIFESTYLE Community Bed	Spare Part
2	2	GSP392212	Recline Linkages for Backrest Frame Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
<b>3+9+D</b>	1	GSP392196	Backrest Assy - Complete - Incl pos 3+9+D - Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
4	1	#N/A	Not Available	Spare Part
5	1	#N/A	Not Available	Spare Part
6	1	GSP392207	Platform Assembly Backrest Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
7	1	GSP392208	Platform Assembly Foot End / Knee Bend Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
9	1	#N/A	Not Available	Spare Part
12	1	#N/A	Not Available	Spare Part
14	2	#N/A	Not Available	Spare Part
15	2	GSP392340	Bed End Shrouds Sublime Teak Aspire MAXI LIFESTYLE Community Bed	Spare Part
17	1	GSP392208	Platform Assembly Foot End / Knee Bend Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
18	1	GSP392285	Locking Latch for Folding Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
25	8	GSP392232	Common Hardware	Spare Part
26	12	#N/A	Common Hardware	Spare Part
27	1	GSP392282	Locking Pin for Backrest Aspire MAXI LIFESTYLE Community Bed	Spare Part
29	2	#N/A	Not Available	Spare Part
30	2	#N/A	Not Available	Spare Part
31	16	#N/A	Common Hardware	Spare Part
34	4	GSP392232	Slider Bush for Platform Frame Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
35	2	GSP392224	End Cap Kit For Mattress PlatForm Aspire MAXI LIFESTYLE Community Bed	Spare Part
36	2	GSP392230	Sliding Block for Backrest Aspire MAXI LIFESTYLE / LIFESTYLE Community Care Bed	Spare Part
38	2	GSP392275	Rubber Caps White Aspire MAXI LIFESTYLE / LIFESTYLE Community Care Bed	Spare Part
39	16	#N/A	Not Available	Spare Part
42	2	GSP392224	End Cap Kit For Mattress PlatForm Aspire MAXI LIFESTYLE Community Bed	Spare Part
43	2	#N/A	Not Available	Spare Part
45	6	GSP392223	End Cap Kit For Mattress Frame Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
46	4	GSP392223	End Cap Kit For Mattress Frame Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
48	4	#N/A	Sourced Locally / Or Cut a Piece of Foam Pad	Spare Part
49	8	#N/A	Common Hardware	Spare Part
51	4	GSP392326	Castor Total Locking Aspire MAXI LIFESTYLE Community Bed	Spare Part
52	10	#N/A	Common Hardware	Spare Part
53	6	#N/A	Common Hardware	Spare Part
54	4	GSP392232	Slider Bush for Platform Frame Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
55	10	#N/A	Common Hardware	Spare Part
56	2	#N/A	Common Hardware	Spare Part
57	4	#N/A	Common Hardware	Spare Part
58	4	GSP392326	Castor Total Locking Aspire MAXI LIFESTYLE Community Bed	Spare Part
59	8	#N/A	Common Hardware	Spare Part
60	1	GSP392285	Locking Latch for Folding Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
61	4	GSP392326	Castor Total Locking Aspire MAXI LIFESTYLE Community Bed	Spare Part
62	6	#N/A	Common Hardware	Spare Part

63	17	#N/A	Common Hardware	Spare Part
64	1	GSP392307	Sticker Pack Aspire MAXI LIFESTYLE Community Care Bed	Spare Part
65	2	GSP392337	Spring Push Pin for Side Rails Aspire MAXI LIFESTYLE Community Bed	Spare Part
20	1	GSP392145	Control Box Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
20	1	GSP392150	Mains Power Cable Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
21	1	GSP392127	Back Actuator Linak Aspire MAXI LIFESTYLE Community Care Bed	Electronics
22	1	GSP392132	Foot Actuator Linak Aspire MAXI LIFESTYLE Community Care Bed	Electronics
21/22	1	GSP392135	Cable Back / Foot Actuator Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
23	3	GSP392025	Actuator Mounting Hardware (Quick Release Pin) Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
24	2	GSP392122	Head / Foot End Actuator Hi Lo Linak Aspire MAXI LIFESTYLE Community Care Bed	Electronics
24a	1	GSP392140	Cable Foot End Actuator Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
24a	1	GSP392142	Cable Head End Actuator Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
28	3	GSP392025	Actuator Mounting Hardware (Quick Release Pin) Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
32	5	GSP392026	Actuator Mounting Hardware (Quick Release Clip) Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
33	5	GSP392026	Actuator Mounting Hardware (Quick Release Clip) Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
40	1	GSP392155	Handset Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
41	1	GSP392160	Plastic Key for Handset Linak Aspire MAXI LIFESTYLE / LIFESTYLE / COMMUNITY CARE Bed	Electronics
66	2	GSP392195	Folding Bracket Assembly – Aspire MAXI LIFESTYLE COMMUNITY CARE Bed	Spare Parts