


Aspire Lift Recline Chair MONTANA - Spare Parts Diagram



| Product Code | Product Description                                                           |
|--------------|-------------------------------------------------------------------------------|
| CHP227400    | Chair Powerlift Recline - Aspire Montana Maxi - Bronze - Single Motor - 220kg |

| No. | Item Number | Item Description                                                                                    |
|-----|-------------|-----------------------------------------------------------------------------------------------------|
| 1   | GSP510060   | Actuator - Single Motor - Aspire Lift Recline Chair MONTANA                                         |
| 1   | GSP510061   | Actuator - Single Motor - 6000N - Aspire Lift Recline Chair MONTANA - High SWL                      |
| 2   | GSP511250   | Base Frame - Aspire Lift Recline Chair MONTANA - MAXI                                               |
| 3   | GSP510190   | Castor - Standard - Aspire Lift Recline Chair MONTANA                                               |
| 4   | GSP510454   | Cover Set - Armrest / Headrest - Aspire Lift Recline Chair MONTANA                                  |
| 5   | GSP510454   | Cover Set - Armrest / Headrest - Aspire Lift Recline Chair MONTANA                                  |
| 6   | GSP510350   | Plastic Foot - Aspire Lift Recline Chair MAUI/MONTANA/RAPHAEL / RENOIR/VERONA                       |
| 7   | GSP510355   | Plastic End Cap - Side Frame - Aspire Lift Recline Chair MAUI/MONTANA/RAPHAEL / RENOIR/VERONA       |
| 8   | GSP510391   | Handset - Aspire Lift Recline Chair MONTANA/VERONA - NEW - Aspire Print                             |
| 9   | GSP510404   | Mains Power Plug (AU) - Aspire Lift Recline Chair MAUI/MONTANA/VERONA                               |
| 10  | GSP510411   | Power Pack / Transformer - Aspire Lift Recline Chair MONTANA/VERONA/MAUI                            |
| 11  | GSP511255   | Seat Extension Frame - Aspire Lift Recline Chair MONTANA - MAXI                                     |
| 12  | GSP400648   | Transformer ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds |
| 13  | GSP400649   | Battery ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds     |
| 14  | GSP510442   | Seat Detachable - Aspire Lift Recline Chair MONTANA                                                 |
| 15  | GSP511350   | Backrest Frame - Aspire Lift Recline Chair MONTANA - MAXI                                           |




12/13



3



2



6



7



9



10



8




4



5



15



11



14



1