


Aspire Powerlift Recline - IDAHO / OREGON - Spare Parts Diagram



Product Code	Product Description
CHP227150	Chair Powerlift Recline - Aspire OREGON - Charcoal Fabric - Space Saver - 150kg ^
CHP227155	Chair Powerlift Recline - Aspire OREGON - Latte Vinyl - Space Saver - 150kg ^
CHP227160	Chair Powerlift Recline - Aspire IDAHO - Charcoal Fabric - Dual Action - 150kg ^
CHP227165	Chair Powerlift Recline - Aspire IDAHO - Latte Vinyl - Dual Action - 150kg ^

Spare Parts		
No.	Item Number	Item Description
1	GSP511105	Actuator - Backrest - Aspire IDAHO Recline Chair
2	GSP511100	Actuator - Lift / Legrest - Aspire IDAHO Recline Chair
2	GSP511107	Actuator - Lift / Recline - Aspire OREGON Recline Chair
3	GSP511103	Actuator Mounting Pins - Aspire IDAHO / OREGON Recline Chair
4	GSP511155	Backrest Brackets - Aspire IDAHO / OREGON Recline Chair - PAIR
5	GSP511150	Castor - Locking - Aspire IDAHO / OREGON Recline Chair
5A	GSP511152	Castor - Non Locking - Aspire IDAHO / OREGON Recline Chair
5B	GSP511160	Castor Bracket - w/mounting hardware - Aspire IDAHO / OREGON Recline Chair
6	GSP511110	Handset - Aspire IDAHO Recline Chair
7	GSP511112	Handset - Aspire OREGON Recline Chair
8	GSP511115	Mains Power Cable / Transformer - Aspire IDAHO/ OREGON Recline Chair
9	GSP511125	Plastic End Cap - Side Frame - Aspire IDAHO / OREGON Recline Chair
10	GSP511120	Plastic Foot - Aspire IDAHO / OREGON Recline Chair
14	GSP400649	Battery ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds
14	GSP400648	Transformer ONLY - for Battery Backup - All Aspire Lift Recline Chair / ComfiMotion Care / LUX Beds

Fabric Covers		
No.	Item Number	Item Description
11/12	GSP511130	Cover - Armrest / Headrest - Aspire IDAHO / OREGON Recline Chair - Latte Vinyl
11/12	GSP511140	Cover - Armrest / Headrest - Aspire IDAHO / OREGON Recline Chair - Charcoal Fabric
13	GSP511220	Base Frame - Aspire IDAHO Recline Chair
13	GSP511310	Base Frame - Aspire OREGON Recline Chair