


Spare Parts Breakdown - Aspire Over Toilet Aid



BTT145000	Over Toilet Aid - Aspire - 460mm - Treated Steel - 175kg
BTT145005	Over Toilet Aid - Aspire - 460mm - Aluminium - 175kg
BTT145010	Over Toilet Aid - Aspire - Wide - 525mm - Treated Steel - 175kg
BTT145015	Over Toilet Aid - Aspire - Wide - 525mm - Aluminium - 175kg

1	GSP390800	Rubber Ferrule/Tip
2	BTT145100	Clip on Toilet Seat
2	BTT145110	Clip on Padded Toilet Seat (not pictured)
3	GSP390802	Height Adjustable Leg (Single Holes)
4	GSP390803	Height Adjustable Leg (Double Holes)
5	GSP390797	Push Button & Spring
6	GSP390796	Plastic Leg Insert
7	GSP390801	Plastic Silencer
8	BTT145105	Splash Guard