

## Spare Parts Manual

---




Item Number	Item Description
BTS148112	Aspire Shower Trolley - 200kg

**This manual is for internal use only**

**1 Electronics**

Ref	Item Number	Item Description	Image	Type




Ref	Item Number	Item Description	Image	Type
1	GSP146704	Control Box - Aspire Shower Trolley		Spare Part
2	GSP146700	Actuator - HILO - Aspire Shower Trolley		Spare Part
3	GSP146702	Actuator - Backrest / Tilt - Aspire Shower Trolley		Spare Part
4	GSP146710	Handset - Aspire Shower Trolley		Spare Part

1 Electronics				
Ref	Item Number	Item Description	Image	Type
5	MBT401050	Battery - Aspire Shower Trolley		Spare Part
6	MBT401055	Charging Cradle - w/Transformer Cable - Aspire Nova Standing / Aura Floor Lifter / Shower Trolley		Spare Part
6A	MBT401060	Power Transformer - w/Cable - for Charging Cradle - Aspire Nova Standing / Aura Floor Lifter / Shower Trolley		Spare Part

**2 Base Frame**

Ref	Item Number	Item Description	Image	Type






Ref	Item Number	Item Description	Image	Type
7	GSP146712	Castor - Total Lock - Aspire Shower Trolley		Spare Part
8	GSP146714	Castor - Directional Lock Aspire Shower Trolley		Spare Part
9	GSP146720	Base Frame - Complete - w/Castors - Aspire Shower Trolley		Spare Part
10	GSP146722	Stopper Block - Base Frame - Aspire Shower Trolley		Spare Part

2 Base Frame				
Ref	Item Number	Item Description	Image	Type
11	GSP146724	Nylon Roller - w/mounting Hardware - on Scissor Foot End - Aspire Shower Trolley		Spare Part
12	GSP146726	HILO Scissor - Complete - w/mounting hardware and Rollers - Aspire Shower Trolley		Spare Part
13	GSP146728	Hardware Set - for Scissor Hinge Points - Aspire Shower Trolley		Spare Part

**3 Deck and under frame**

Ref	Item Number	Item Description	Image	Type



Ref	Item Number	Item Description	Image	Type
14	GSP146734	Deck Support Frame - Aspire Shower Trolley		Spare Part
15	GSP146736	Stopper Block - w/mounting hardware - Under Deck - Aspire Shower Trolley		Spare Part
17	GSP146738	Support Roller - w/mounting hardware - Backrest Deck - Aspire Shower Trolley		Spare Part
18	GSP146740	Backrest Deck Assembly - w/Mounting Hardware - Aspire Shower Trolley		Spare Part

3 Deck and under frame				
Ref	Item Number	Item Description	Image	Type
20	GSP146742	Main Deck Assembly - w/Mounting Hardware - Aspire Shower Trolley		Spare Part
22	BTS148140	Extension KIT - Complete - w/Hardware - Aspire Shower Trolley		Accessory
28	GSP146746	Foot End Link - w/Hardware - Pair - Aspire Shower Trolley		Spare Part
29	GSP146748	Head End Link - w/Hardware - Pair Aspire Shower Trolley		Spare Part
30	GSP146750	Twist Clip - to hold Liner Aspire Shower Trolley		Spare Part
32	GSP146752	Backrest / Tilt Bar - w/Head Connection & mounting hardware - Aspire Shower Trolley		Spare Part

**4 Siderails & End Rails**


Ref	Item Number	Item Description	Image	Type

Ref	Item Number	Item Description	Image	Type
33	GSP146760	Siderail - Plastic Rail only - (1 pce) - LHS/RHS reversible - Aspire Shower Trolley		Spare Part
34	GSP146764	Frame Assly - Side Rail - Complete - (1 pce) - LHS/RHS reversible - Aspire Shower Trolley		Spare Part
35	GSP146762	Head / Footend Rail - Plastic Rail only - (1 pce) - Head/Foot reversible - Aspire Shower Trolley		Spare Part
36	GSP146766	Frame Assly - Head/Footend Rail - (1 pce) - Head/Foot reversible - Aspire Shower Trolley		Spare Part

4 Siderails & End Rails				
Ref	Item Number	Item Description	Image	Type
40	GSP146768	Spring - for Siderail Hinge - Aspire Shower Trolley		Spare Part
41	GSP146770	Siderail Release Bar - Aspire Shower Trolley		Spare Part

## 6 Liner, Drain and Accessories

Ref	Item Number	Item Description	Image	Type
Ref	Item Number	Item Description	Image	Type
48	GSP146776	Liner - Standard - w/Pillow, Fittings and Drain Hose - Aspire Shower Trolley		Spare Part
48	BTS148120	Liner - Long - w/Pillow, Fittings and Drainage Hose - Aspire Shower Trolley		Accessory
49	GSP146782	Pillow - Aspire Shower Trolley		Spare Part
50	GSP146778	Drain Coupling & Hose - Aspire Shower Trolley		Spare Part
52	BTS148125	Shower Head Holder - Aspire Shower Trolley		Accessory

53	BTS148130	Accessory Basket - Aspire Shower Trolley	 An illustration showing a grey plastic accessory basket with a handle and a latch, mounted on the side of a white shower trolley. The basket is designed to hold small items like soap or shampoo.	Accessory
----	-----------	--	---	-----------