

Aspire VIDA AP (formerly DASH) Wheelchairs - Spare Parts Manual



Product Code	Product Description
MWS449850	WChair Folding - Aspire DASH AP - 450mm - 135kg - Silver - Discontinued
MWS449855	WChair Folding - Aspire DASH AP - 450mm - 135kg - Other Colour - Discontinued
MWS449874	WChair Folding - Aspire VIDA - AP - 400mm - 135kg - Mustard
MWS449878	WChair Folding - Aspire VIDA - AP - 450mm - 135kg - Teal
MWS449882	WChair Folding - Aspire VIDA - AP - 450mm - 135kg - Magenta
MWS449886	WChair Folding - Aspire VIDA - AP - 500mm - 135kg - Olive

Old Dash W/Chair Frames are Silver and New VIDA W/Chair Frames are Gunmetal Grey

No.	Item Number	Item Description
1,2	GSP489105	Backrest Cane with Handgrip - Aspire DASH AP / SP Wheelchair - PAIR
1,2	GSP489146	Backrest Cane w/Handgrip - New G/Grey Frame - Aspire VIDA AP/SP WChair - PAIR
3	GSP489100	Armrest Pad - Aspire DASH AP / SP Wheelchair - PAIR
4	GSP489125	Mounting Hardware - for Backrest Cane - Aspire DASH AP / SP Wheelchair - PAIR
5	N/A	NOT AVAILABLE
6	GSP489150	Side Arm Fabric Guard - Aspire DASH AP / SP Wheelchair - PAIR
7	GSP489130	Rear Wheel Assembly - Aspire DASH AP Wheelchair
7A	GSP489172	Fixed Axle - (Nut, Bolt and Washer) - Aspire DASH AP / SP Wheelchair
8	GSP489110	Brake Handle and Cable Assembly - Aspire DASH AP Wheelchair - PAIR
9	GSP489120	Front Fork and Castor Assembly - Aspire DASH AP / SP Wheelchair - LEFT
9	GSP489121	Front Fork and Castor Assembly - Aspire DASH AP / SP Wheelchair - RIGHT
10	GSP489155	S/A Footrest Assembly - Aspire DASH AP / SP Wheelchair - PAIR
10	GSP489148	S/A Footrest Assembly - New G/Grey Frame - Aspire VIDA AP/SP WChair - PAIR
10A	GSP489160	Foot Plates and Calf Strap - Aspire DASH AP / SP Wheelchair - LEFT
10A	GSP489162	Foot Plates and Calf Strap - Aspire DASH AP / SP Wheelchair - RIGHT
11	GSP489140	Seat and Back Upholstery - Aspire DASH AP / SP Wheelchair - Green
11	GSP489145	Seat and Back Upholstery - Aspire DASH AP / SP Wheelchair - Orange
11	GSP489141	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 400mm - Mustard
11	GSP489142	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 450mm - Teal
11	GSP489143	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 450mm - Magenta
11	GSP489144	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 500mm - Olive
12	N/A	NOT AVAILABLE
13	GSP489165	End Cap - Aspire DASH AP / SP Wheelchair
14	WCA848275	WChair Anti Tip Bar - Aspire VIDA (Formerly DASH) Wheelchair

SP Wheelchairs Next page

Aspire VIDA SP (formerly DASH) Wheelchairs - Spare Parts Manual



Product Code	Product Description
MWS449860	WChair Folding - Aspire DASH SP - 450mm - 135kg - Silver - Discontinued
MWS449865	WChair Folding - Aspire DASH SP - 450mm - 135kg - Other Colour - Discontinued
MWS449876	WChair Folding - Aspire VIDA - SP - 400mm - 135kg – Mustard
MWS449880	WChair Folding - Aspire VIDA - SP - 450mm - 135kg – Teal
MWS449884	WChair Folding - Aspire VIDA - SP - 450mm - 135kg - Magenta
MWS449888	WChair Folding - Aspire VIDA - SP - 500mm - 135kg - Olive

Old Dash W/Chair Frames are Silver and New VIDA W/Chair Frames are Gunmetal Grey

No.	Item Number	Item Description
1,2	GSP489105	Backrest Cane with Handgrip - Aspire DASH AP / SP Wheelchair - PAIR
1,2	GSP489146	Backrest Cane w/Handgrip - New G/Grey Frame - Aspire VIDA AP/SP WChair - PAIR
3	GSP489100	Armrest Pad - Aspire DASH AP / SP Wheelchair - PAIR
4	GSP489125	Mounting Hardware - for Backrest Cane - Aspire DASH AP / SP Wheelchair - PAIR
5	N/A	NOT AVAILABLE
6	GSP489150	Side Arm Fabric Guard - Aspire DASH AP / SP Wheelchair - PAIR
7	GSP489135	Rear Wheel Assembly - Aspire DASH SP Wheelchair
7A	GSP489170	Quick Release Axle - Aspire DASH SP Wheelchair
8	GSP489115	Brake Handle and Cable Assembly - Aspire DASH SP Wheelchair - PAIR
9	GSP489120	Front Fork and Castor Assembly - Aspire DASH AP / SP Wheelchair - LEFT
9	GSP489121	Front Fork and Castor Assembly - Aspire DASH AP / SP Wheelchair - RIGHT
10	GSP489155	Swing Away Footrest Assembly - Aspire DASH AP / SP Wheelchair - PAIR
10	GSP489148	S/A Footrest Assembly - New G/Grey Frame - Aspire VIDA AP/SP WChair - PAIR
10A	GSP489160	Foot Plates and Calf Strap - Aspire DASH AP / SP Wheelchair - LEFT
10A	GSP489162	Foot Plates and Calf Strap - Aspire DASH AP / SP Wheelchair - RIGHT
11	GSP489140	Seat and Back Upholstery - Aspire DASH AP / SP Wheelchair - Green
11	GSP489145	Seat and Back Upholstery - Aspire DASH AP / SP Wheelchair - Orange
11	GSP489141	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 400mm - Mustard
11	GSP489142	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 450mm - Teal
11	GSP489143	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 450mm - Magenta
11	GSP489144	Seat and Back Upholstery - Aspire VIDA AP/SP Wheelchair - 500mm - Olive
12	N/A	NOT AVAILABLE
13	GSP489165	End Cap - Aspire DASH AP / SP Wheelchair
14	WCA848275	WChair Anti Tip Bar - Aspire VIDA (Formerly DASH) Wheelchair