


Spare Parts Diagram
 BEB050140 / BEB050150
 Hospital Bed - Malvestio Delta 4 Ward Bed - 300kg
 Hospital Bed - Malvestio Delta 4 Ward Bed W/ Elec Drive Wheel & Weigh Scale - 300kg

* For exploded diagram and a more comprehended list please check the service manuals located on the intranet which can be found by clicking the link here:
<http://intranet.aidacare.com.au/sites/rm/Spare%20Parts%20Diagrams/Malvestio%20Hospital%20Equipments/Delta%204%20Bed/Malvestio%20Delta%204%20Bed%20-%20Service%20Manual%20-%2020210421.pdf>

Aidacare Code	Reference	Supplier Code	Aidacare Description
GSP497300	Page 61, H	871007	Actuator - Backrest - TiMotion - Malvestio Delta 4 Ward Bed
GSP497610		872208	Cable - Backrest Actuator - TiMotion - Malvestio Delta 4 Ward Bed
GSP497301	Page 61, H	871027	Actuator - Backrest - Linak - Malvestio Delta 4 Ward Bed
GSP497611		872251	Cable - Backrest Actuator - Linak - Malvestio Delta 4 Ward Bed
GSP497305	Page 61, I	871006	Actuator - Knee Break - TiMotion - Malvestio Delta 4 Ward Bed
GSP497615		872209	Cable - Knee Break Actuator - TiMotion - Malvestio Delta 4 Ward Bed
GSP497306	Page 61, I	871028	Actuator - Knee Break - Linak - Malvestio Delta 4 Ward Bed
GSP497616		872252	Cable - Knee Break Actuator - Linak - Malvestio Delta 4 Ward Bed
GSP497310	Page 61, F/G	871005	Actuator - Hi/Lo - TiMotion - Malvestio Delta 4 Ward Bed
GSP497600		872206	Cable - Hi/Lo Actuator - Head end - TiMotion - Malvestio Delta 4 Ward Bed
GSP497605		872207	Cable - Hi/Lo Actuator - Foot end - TiMotion - Malvestio Delta 4 Ward Bed
GSP497311	Page 61, F/G	871026	Actuator - Hi/Lo - Linak - Malvestio Delta 4 Ward Bed
GSP497601		872248	Cable - Hi/Lo Actuator - Head end - Linak - Malvestio Delta 4 Ward Bed
GSP497606		872249	Cable - Hi/Lo Actuator - Foot end - Linak - Malvestio Delta 4 Ward Bed
GSP497315	NA	RC99251	Battery Backup - Twin 1.2Amps - Malvestio Delta 4 Ward Bed
GSP497320	Page 67, 28	41000029	Bed End Support Frame - Malvestio Delta 4 Ward Bed
GSP497325	Page 67, 22	41000028	Bed Extension Frame - Malvestio Delta 4 Ward Bed
GSP497328	Page 68, 69	41000024	Backrest Frame - Malvestio Delta 4 Ward Bed
GSP497330	Page 73, 2	826642	Side Rail - Left Foot - Malvestio Delta 4 Ward Bed (Plastic Shell Only)
GSP497335	Page 73, 1	826703	Side Rail - Left Head - Malvestio Delta 4 Ward Bed (Plastic Shell Only)
GSP497340	Page 73, 3	826574	Side Rail - Right Foot - Malvestio Delta 4 Ward Bed (Plastic Shell Only)
GSP497345	Page 73, 4	826702	Side Rail - Right Head - Malvestio Delta 4 Ward Bed (Plastic Shell Only)
GSP497350	Page 62, 5	825832	Brake Cam - Castor - Malvestio Delta 4 Ward Bed
GSP497355	Page 62, 3	826713	Brake Pedal - Malvestio Delta 4 Ward Bed
GSP497360	Page 65, 11	825698	Stopper - Frame - Malvestio Delta 4 Ward Bed
GSP497365	Page 60, A	497843_230AUS	Control Box and Mains Power Cable - Malvestio Delta 4 Ward Bed - Malvestio Delta 4 Ward Bed
GSP497100	Page 61, E	870647	Control Panel - Nurse Attendant - Malvestio Delta 4 Ward Bed
GSP497375	Page 61, D	880620	Control Panel - Side Rail - Left - Malvestio Delta 4 Ward Bed
GSP497380	Page 61, C	880613	Control Panel - Side Rail - Right - Malvestio Delta 4 Ward Bed
GSP497390	Page 69, 1	826550	Mattress Platform Cover - Backrest - Malvestio Delta 4 Ward Bed
GSP497395	Page 62, 2	826552	Shroud Cover - Base Frame Foot - Malvestio Delta 4 Ward Bed
GSP497400	Page 62, 1	826553	Shroud Cover - Base Frame Head - Malvestio Delta 4 Ward Bed
GSP497405	Page 69, 2/3/6	826594	Mattress Platform Cover - Central Section - Malvestio Delta 4 Ward Bed
GSP497410	Page 2, 2 / 3	826560/825536	Corner Bumper and Insert Assembly - Corner Frame for Head / Foot Board - Malvestio Delta 4 Ward Bed
GSP497415	Page 69, 4	826548	Mattress Platform Cover - Knee Break - Malvestio Delta 4 Ward Bed
GSP497420	Page 69, 5	826547	Mattress Platform Cover - Leg Section - Malvestio Delta 4 Ward Bed

GSP497425	Page 69, 22 / Page 76, 5	879129 / 879145	CPR Handle and Release Cable Assembly - Malvestio Delta 4 Ward Bed
GSP497430	NA	870948	Fuse - 10Amps - Malvestio Delta 4 Ward Bed
GSP497435	Page 61, Q	870977	Handset - Malvestio Delta 4 Ward Bed
GSP497440	Page 65, 1	49000188	Head / Foot Board - Malvestio Delta 4 Ward Bed
GSP872608	NA	49000314	Head Board - for Drive Wheel - AB 430 Blu Iris - Malvestio Delta 4 / Aria Bed
GSP497445	NA	826451	Hook - for Mains Power Cable - Malvestio Delta 4 Ward Bed
GSP497447	NA	880179	Mains Power Cable - for Control Box - Malvestio Delta 4 Ward Bed
GSP497455	Page 61, S / T	872079	Brake Microswitch LED - Malvestio Delta 4 Ward Bed
GSP497460	NA	870358	PCB Board - - Malvestio Delta 4 Ward Bed
GSP497465	Page 69, 23	879134	Quick Release Cable - Backrest - TiMotion - LEFT - Malvestio Delta 4 Ward Bed
GSP497470	Page 69, 23	879132	Quick Release Cable - Backrest - TiMotion - RIGHT - Malvestio Delta 4 Ward Bed
GSP497475	Page 69, 23	879150	Quick Release Cable - Backrest - Linak - RIGHT/LEFT - Malvestio Delta 4 Ward Bed
GSP497480	Page 65, 23	43000004	Release Handle - For Bed Extension - Malvestio Delta 4 Ward Bed
GSP497485	Page 73, 25	872269	Shock Absorber - Rotary - Malvestio Delta 4 Ward Bed
GSP497495	Page 65, 17	826086	Universal Hook - for Bar - Malvestio Delta 4 Ward Bed
MBT410025		MKES12012	Battery - 12 Amp/Hour - MK Sealed Gel - for Power Wheelchairs/Scooters & Malvestio Drive Wheel Bed
GSP497346	1	RC99702	Side Rail - with mounting mechanism - Right Head - Ref - 337110 - Delta 4 Bed
GSP497347	2	RC99703	Side Rail - with mounting mechanism - Left Head - Ref - 337110 - Delta 4 Bed
GSP497348	3	RC99706	Side Rail - with mounting mechanism - Right Foot - Ref - 337110-337160 - Delta 4 Bed
GSP497349	4	RC99707	Side Rail - with mounting mechanism - Left Foot - Ref - 337110-337160 - Delta 4 Bed



GSP497497	Page 72, 73	RC99662	Unhooking Handle Kit for Side Rails (826618 Hook) - for Malvestio Delta 4 / Aria Beds
-----------	-------------	---------	---



BEA825094	Page 48	337027	Folding Monitor Tray - Malvestio Delta 4 and Aria Bed	
				
GSP046310	Page 73	RC99786	Degree Indicator - Inclinometer Kit - Right Head Side - Malvestio Delta 4 / Aria Bed	
GSP046312		RC99787	Degree Indicator - Inclinometer Kit - Left Head Side - Malvestio Delta 4 / Aria Bed	
GSP046314		RC99788	Degree Indicator - Inclinometer Kit - Right Foot Side - Malvestio Delta 4 / Aria Bed	
GSP046316		RC99789	Degree Indicator - Inclinometer Kit - Left Foot Side - Malvestio Delta 4 / Aria Bed	
