

## Malvestio - Idea 2 - Treatment Chair

### Spare Parts Manual

---





Item Number	Item Description
CHT100710	Treatment Chair - Malvestio IDEA 2 - Navy Blue Vinyl - 200kg


**This manual is for internal use only**

**1 Spare Parts**




Ref	Item Number	Item Description	Image	Type
1	GSP872696	Twin Wheel - 125mm - Braking - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
1	GSP872698	Twin Wheel - 125mm - Antistatic Directional Lock - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
2	GSP794710	Padded Backrest - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
3	GSP794712	Upholestered Seat - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part






1 Spare Parts				
4	GSP794714	Padded Legrest - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
5	GSP794706	Padded Headrest - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
6	GSP794702	Armrest Assembly - <b>Pair</b> - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
7	GSP794708	Armrest Pad - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
8	GSP794700	Plastic Brake Pedal - for Wheels - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
9	GSP794716	Flap Closure - (Backrest/Seat) - Zipper - Malvestio IDEA 2 Treatment Chair (384420)	N/A	Spare Part

1 Spare Parts				
9A	GSP794704	Backrest Insert - Malvestio IDEA 2 Treatment Chair (384420)		Spare Part

**2 Electronics**



Ref	Item Number	Item Description	Image	Type
10	GSP872660	Power Cable w/Australian Plug - for Malvestio Treatment Chairs (384365_AU / 384420)	Not on Drawing	Spare Part
11	GSP872680	Control Box - 497863 CL.1 230V - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
12	GSP872682	Checkup CPU for Control Box - 4M CLASS I - for Malvestio IDEA 2 Treatment Chair (384420) <i>(this is a CPU inside the Control Box)</i>		Spare Part
13	GSP872684	Lifting Column - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part

2 Electronics				
14	GSP872686	Handset - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
15	GSP872688	Actuator - Backrest - TIMOTION 870663 - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
16	GSP872690	Actuator - Seat / Trend - TIMOTION 870661 - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
17	GSP872692	Actuator - Legrest - TIMOTION 870662 - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part
18	GSP872694	Battery - Dewert - IP66 + CABLE 0.42M - for Malvestio IDEA 2 Treatment Chair (384420)		Spare Part