

Malvestio - Idea Clinic - Treatment Chair Spare Parts Manual



Item	Description
CHT100700	Treatment Chair - Malvestio IDEA CLINIC - Navy Blue Vinyl - 200kg

This manual is for internal use only

1 Spare Parts




Ref	Item Number	Item Description	Image	Type
A	GSP872662	Castor - Tente - 100mm - Braking - w/Pintle - Malvestio Idea Clinic Treatment Chair		Spare Part
1	GSP872664	Castor - Tente - 100mm - Directional - w/Pintle - Malvestio Idea Clinic Treatment Chair		Spare Part
B	GSP794754	Upholstered Backrest - Malvestio Idea Clinic Treatment Chair (384365_AU)		Spare Part
C	GSP794752	Upholstered Seat - Malvestio Idea Clinic Treatment Chair (384365_AU)		Spare Part




1 Spare Parts				
D	GSP794758	Padded Legrest - Malvestio Idea Clinic Treatment Chair (384365_AU)		Spare Part
E	GSP794756	Upholstered Footrest - Malvestio Idea Clinic Treatment Chair (384365_AU)		Spare Part
F	GSP794760	Padded Headrest - Pillow - Malvestio Idea Clinic Treatment Chair (384365_AU)	N/A	Spare Part
G	GSP794750	Armrest - suits both LHS/RHS - Malvestio Idea Clinic - Treatment Chair		Spare Part
H	GSP794762	Lever - Armrest - GS Unlocking - Malvestio Idea Clinic Treatment Chair (384365_AU)		Spare Part
I	GSP794764	Armrest Gas Strut - 300N - Malvestio Idea Clinic Treatment Chair (384365_AU)	N/A	Spare Part


1 Spare Parts

J	GSP872670	Locking Lever - Armrest Height - Malvestio Idea Clinic Chair		Spare Part
---	-----------	--	--	------------

1 Electronics



Ref	Item Number	Item Description	Image	Type
7	GSP872650	Handset - 4 Function / 8 Button - Malvestio Idea Clinic - Treatment Chair		Spare Part
8	GSP872656	Control Box - Malvestio Idea Clinic - Treatment Chair		Spare Part
9	GSP872654	Backrest Actuator - Malvestio Idea Clinic - Treatment Chair		Spare Part
10	GSP872652	Seat Actuator - Malvestio Idea Clinic - Treatment Chair		Spare Part

1 Electronics				
11	GSP872658	Lifting Column - Malvestio Idea Clinic - Treatment Chair		Spare Part
12	GSP872660	Power Cable w/Australian Plug - for Malvestio Treatment Chairs (384365_AU / 384420)	Not on Drawing	Spare Part
13	GSP872668	Battery back up - TiMotion - for Malvestio Idea Clinic Treatment Chair	Not on Drawing	Spare Part