

Mangar Rhino/Swift and Airflo Duo Pump Spare Parts Manual



Product code: CE1000
Airflo Duo




Product Code: MPCA241325
Rhino & Swift

Item Number	Item Description
GSP874504	Replacement Pump - Airflo Duo - Mangar RHINO & SWIFT
LSS403830	Lift Cushion - Mangar RHINO & SWIFT Supine System - 500kg

This manual is for internal use only

1 Airflo Duo Pump



Ref	Item Number	Supplier Code	Description	Type
3	GSP874505	AT0010AU	Mains Power Adaptor - Airflo Duo - Mangar RHINO & SWIFT	Spare Part
7	GSP874504	CE1023	Replacement Pump - Airflo Duo - Mangar RHINO & SWIFT	Spare Part
9	GSP874506	AT0014AU	Battery Charger - Airflo Duo - Mangar RHINO & SWIFT	Spare Part
10	GSP874507	AT0015	Hose Kit - Mangar - Airflo Duo - Mangar RHINO & SWIFT	Spare Part
11	GSP874501	CE1024	Replacement Battery Pack - Airflo Duo - Mangar RHINO & SWIFT	Spare Part
12	POA	CE1025	Airflo Duo Filter Kit	Spare Part
13	POA	CE1027	Airflo Duo Motor Kit	Spare Part
14	POA	CE1028	Airflo Duo PCB Kit	Spare Part
15	POA	CE1029	Airflo Duo Front Case Kit	Spare Part
16	POA	CE1030	Airflo Duo Rear Case Kit	Spare Part

1 Airflo Duo Pump				
17	POA	AT0006	Airflo Duo -Mains Power Supply	Spare Part
18	POA	CE1057	Airflo Duo Rear Hooks Kit	Spare Part

2 Rhino & Swift Cushions



30



31



3



33



34



HJ1114
Valve Removal Tool Kit

35



36

Ref	Item Number	Supplier Code	Description	Type
30	GSP874518	HJ1025	STS Restraint Strap Kit - Mangar RHINO & SWIFT	Spare Part
31	GSP874519	HG0008	Rhino Handle Kit - Mangar RHINO & SWIFT	Spare Part
32	GSP874520	HG1024	Swift Side Pull Strap Kit - Mangar RHINO & SWIFT	Spare Part
33	POA	HJ1113	Rhino Valve and Dust Cap Kit - Rhino	Spare Part
34	POA	HJ1112	Swift Valve Kit - Swift	Spare Part
35	POA	HJ1114	Valve Removal Tool Kit - Rhino & Swift	Spare Part
36	POA	AT0012	Supine Transfer System Bag Kit - Rhino	Spare Part
36	POA	AT0016	Swift Bag Kit - Swift	Spare Part