

MASS QLD ORDER FORM SUMMARY PAGE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

DEALER INFORMATION

Date	_____	Name	_____
Account Code	_____	Address	_____
Sales Rep	_____	City / State	_____
Mobile	_____	Postcode	_____
Signature	_____	Fax	_____
		Email	_____

CUSTOMER INFORMATION

Customer Name	_____	Contact	_____
MASS Plaque #	_____	Address	_____
MASS PO#	_____	City / State	_____
Signature	_____	Postcode	_____
Rep Name & Ph #	_____	Phone	_____
		Fax	_____
		Email	_____

PRESCRIBER INFORMATION

Prescriber Name	_____	Organisation	_____
Mobile	_____	Address	_____
		City / State	_____
		Postcode	_____
Signature	_____	Phone	_____
		Fax	_____
		Email	_____

PRODUCT INFORMATION

The OSD Aspire EVOKE is a highly adjustable lightweight aluminum wheelchair designed to provide a personalised fit and feel for every user. Extensive seat and wheel adjustability offers maximum independence for the user without compromising clinical outcomes. With a comprehensive range of accessories available, the EVOKE provides maximum performance and is suitable for indoor and outdoor use.

Max. User Weight	S.W.L 130kg
	15kg (chair weight)
Weight	8.5kg (transport weight)
Seat Height - adjustable (floor to top of fabric)	FRONT: 50cm 47cm 45cm REAR: 45cm 43cm 41cm
Seat Width	35cm (14") 40cm (16") 45cm (18") 50cm (20")
Seat Depth - adjustable	40cm (16") 43cm (17") 45cm (18") 48cm (19") 50cm (20") 53cm (21")
Backrest Height - adjustable (top of seat to top of backrest handle)	40cm - 45cm
Overall Width	Seat width plus 20cm
Armrest Height (top of seat to top of armrest pad)	22cm fixed type, variable height and depth on adjustable option
Approvals	TGA# ARTG224574 AS/NZS tested and approved



MASS QLD ORDER FORM

Effective:

ORDERING INFORMATION: Phone 1300 133 120

ACCOUNT INFORMATION

Invoice to _____	Deliver to _____
Account # _____	Name _____
Address _____	Address _____
_____	_____
Email _____	Email _____
Phone _____	Phone _____
Order Date _____	Marked for _____
PO # _____	Order # _____
Number Required _____	_____

EVOKE STANDARD MASS PACKAGE

Colour	Silver	Seat	Standard padded nylon	Footplate	Angle adjustable footplate
Backrest	Tension adjustable	Back Upholstery	Tension adjustable back	Castor	8"
Rear Wheels	24"	Armrest	22cm	Brake	Push to lock
		Legrest	70° swing in / out	Accessories	Tool kit

FOR OPTIONAL EXTRAS, PLEASE CHOOSE ONE OF THE PACKAGES BELOW:

CHAIR CUSTOMISATION without cushion - PLEASE SELECT

All customisations will be quoted as labour in addition to the base price, customisations can be carried out by Aidacare or by an authorised Aidacare Network Partner.

Seat Width	35 cm	40 cm	45 cm	50 cm		
Seat Depth	40 cm	43 cm	45 cm	48 cm	50 cm	53 cm
Seat Height	<small>FRONT</small> 50 cm	47 cm	45 cm	<small>BACK</small> 45 cm	43 cm	41 cm

ACCESSORIES - PLEASE SPECIFY QUANTITY

Code	Accessory	Price	Qty
WCA858370	Elevating Legrest - Black Left	Not on MASS SOA	\$104
WCA858390	Elevating Legrest - Black Right	Not on MASS SOA	\$104
WCA858440	Stump Support - Black Left		\$45
WCA858450	Stump Support - Black Right		\$45
MSP421020	Armrests (Height & Depth Adj.) - Pair	Not on MASS SOA	\$99
WCA858460	Anti Tip Bars EVOKE 350		\$15
WCA858463	Transport Tie Downs	Not on MASS SOA	\$52
MWT469865	24" Pneumatic Rear Wheels (Pair)	Not on MASS SOA	\$115
WCA858453	Swing Away Legrest EVOKE 350		\$60
WCA858454	Swing Away Legrest EVOKE 400		\$60
WCA858455	Swing Away Legrest EVOKE 450 & 500		\$60

MASS Plaque No. _____

MASS Funded _____

Client Contribution _____

TOTAL PRICE _____