

Quick and Easy Disassembly




Front LED lights shown above
(also includes rear LED light)

With its compact size, the **Sonic** is the perfect choice for your on-the-go lifestyle. This portable and durable scooter quickly dismantles into five lightweight pieces for easy transport and storage, no tools required! With puncture-proof tyres and a detachable battery pack, the Sonic comes loaded with luxury features.



Specifications

Weight Capacity	136 kg (299.8 lbs.)
Turning Radius	110 cm (43.3")
Maximum Speed	Up to 6.8 km/h (4.2 mph)
Overall Length	104 cm (40.9")
Overall Width	49.5 cm (19.4")
Overall Height	92.5 cm (36.4")
Total Weight	45 kg (99 lbs.)
Maximum Range	Up to 20 km (12.4 miles)
Max. Degree of Climb	12°
Motor	24V, 180W
Battery Capacity	20 AH
Battery Weight	6.9 kg (15.2 lbs.) (20 AH)
Charger	2A, 24V
Wheel Size	Front: 20.32 cm (8") Rear: 22.86 (9") (solids)
Ground Clearance	6.5 cm (2.5")
Colour	

Features

- Front and rear LED lights
- Easy access battery compartment
- Adjustable angle tiller with ergonomic delta handle
- Optional rear-view mirror



Rear LED tail lights

NOTE: All specifications subject to change without notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661
(Tel) 03 8770 9600 • (Fax) 03 9703 2720
www.pridemobility.com.au

Pride®
Mobility Products
Australia Pty. Ltd.