








# PRIDE MOBILITY SCOOTERS - CHOICE ASSISTANCE GUIDE

							
<b>Model</b>	<b>Boot Scooter</b>	<b>Go Go LX + CTS</b>	<b>Sonic Scooter</b>	<b>Quest Folding Compact</b>	<b>GT Scooter</b>	<b>Pathrider 10</b>	<b>Pathrider 10 DX</b>
<b>Size</b>	<b>*SMALL*</b>	<b>*SMALL*</b>	<b>*SMALL*</b>	<b>*SMALL*</b>	<b>**MEDIUM**</b>	<b>**MEDIUM**</b>	<b>**MEDIUM**</b>
<b>Code</b>	SCT670613	SCT670612	SCT670665	SCT670630	SCT670670	SCT670615	SCT670620
<b>Range (Max.)</b>	11.5km	17km	20km	20km	30km	40km	33km
<b>Battery</b>	2x 12V 12AH	2x 12V 18AH	2x 12V 20AH	1x 24V 11.5AH Li-Polymer	2 x 12V 35AH	2x 12V 31AH (40AH optional)	2x 12V 31AH (40AH optional)
<b>Seat Type</b>	Low back	Low back	Low Back	Low back	Low Back	Low back	High back
<b>Max Weight</b>	118kg	136kg	136kg	115kg	159kg	159kg	159kg
<b>Std Body Colours</b>	Red	Red	Red	White	Red	Viper Blue	Viper Blue
<b>RRP\$</b>	\$1,455	\$1,981	\$1,299	\$3,052	\$1,999	\$2,512	\$2,998
<b>Features</b>	Lightweight, <b>disassembles</b> , dual-voltage charger, adjustable tiller, flat free tyres, front basket	Lightweight, <b>disassembles</b> , comfort track suspension (CTS), adjustable tiller, front basket	Lightweight, <b>disassembles</b> , LED light, adjustable tiller and delta handle	Lightweight <b>Folding</b> design, folds in three easy steps, adjustable tiller	Advanced suspension, LED light, deluxe tiller, directional signals	High performance, sporty, flat free tyres, back lit battery gauge, armrests, front basket	High performance, sporty, flat free tyres, head & armrests, front & rear suspension

					
<b>Model</b>	<b>Victory 10 LX</b>	<b>Pursuit 130 XL</b>	<b>Pathrider 140 XL</b>	<b>Celebrity Deluxe XL</b>	<b>Sportrider Scooter</b>
<b>Size</b>	<b>**MEDIUM**</b>	<b>***LARGE***</b>	<b>***LARGE***</b>	<b>***LARGE***</b>	<b>***LARGE***</b>
<b>Code</b>	SCT670640	SCT670571	SCT670576	SCT670540	SCT670585
<b>Range (Max.)</b>	26km	25km	35km	40km	50km
<b>Battery</b>	2x 12V 40AH	2x 12V 55AH	2x 12V 70AH	2x 12V 40AH	2x 12V 70AH
<b>Seat Type</b>	High back	CRS Deluxe	Deluxe Captains Seat	Deluxe reclining high back	High back, adjustable
<b>Max Weight</b>	181kg	159kg	181kg	181kg	181kg
<b>Std Body Colours</b>	Candy Apple Red	Candy Apple Red	Candy Apple Red, Silver & Yellow	Candy Apple Red & Silver	Silver
<b>RRP\$</b>	\$3,426	\$3,951	\$5,190	\$5,083	\$6,957
<b>Features</b>	High performance, outdoor use, head & armrests, lumbar support, rear suspension	High performance, sporty, large pneumatic tyres, head & armrests, full lighting package	Top of range sporty, large pneumatic tyres for outdoor, head & armrests, front & rear suspension	High performance, large tyres for outdoor, head & armrests, lumbar support, rear suspension	Sleek outdoor scooter, digital speedometer display, full suspension system, front/rear lighting

## SCOOTER ACCESSORIES

- A. Rear Scooter Bag
- B. Scooter Flag
- C. Sun Canopy
- D. Walking Frame Holder



All information within this document is correct at the time of publication. This document should be used as a guide only. For further enquiries regarding product specifications, please refer to product brochures. \*Range Per Charge varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition & tyre condition. Battery upgrades (where possible) improve Range Per Charge. \*\*All prices are recommended retail. This varies with promotions. Contact Aidacare for the most up to date pricing.