

MANUAL CHAIRS	ADULT								
	Rigid					Folding			Tilt In Space
	QRi	GPV, GPSA	Q7	5R	7R	QXi	Quickie 2 / Q2 Lite	Xenon / Xenon SA / Xenon Dual	SR45
									
Weight (kg)	10	11	6	5.8	4.5	13.25	8/13	4.9	17
Weight Limit (kg)	120	113	120	120	120	120	120 / 120	110 / 125 / 140	120
Seat Width (mm)	300 - 510	300 - 510	300 - 510	300 - 510	300 - 510	360 - 510	310 - 560, 350 - 500	320 - 500	360 - 510
Seat Depth (mm)	300 - 510	300 - 560	300 - 510	300 - 510	300 - 510	360 - 510	310 - 560, 350 - 500	340 - 500	360 - 560
Floor to Seat Ht (mm)	380 - 510	390 - 580	410 - 530	400 - 530	400 - 530	320 - 500	360 - 560	430 - 560	360 - 520
Features	Entry level rigid	Mid range rigid. Swing away leg option. Highly configurable	Top end rigid with extensive configuration	Highly adjustable including centre of gravity and backrest height and angle	Ultra lightweight, freestyle backrest system, highly adjustable	Entry level folding	Highly configurable scripted folding chair	Active user folding with rigid feel. Swing away legs option	Entry level TIS (rotational)

MANUAL CHAIRS	ADULT / PAED	PAEDIATRIC				
	Tilt In Space	Rigid	Folding		Tilt In Space	
	Quickie / Zippie Iris	Zippie Zone	Zippie "GS"	Youngster 3	Zippie X'Cape	Zippie TS
						
Weight (kg)	10/17	6.5	11.5	8.5	8.6	14
Weight Limit (kg)	75 / 158	75	75	85	75	75
Seat Width (mm)	113 / 135 / 158	200 - 460	250 - 460	230 - 410	203 - 406	250 - 460
Seat Depth (mm)	250 - 460, 360 - 560	200 - 480	250 - 510	240 - 410	203 - 406	330 - 560
Floor to Seat Ht (mm)	320 - 500	360 - 480	360 - 510	250 - 410	330 - 482	390 - 500
Features	Highly configurable TIS (rotational)	Active user	Growth system. Front or rear wheel drive	Active user	Xlock Crossbrace for performance, Z-Finity system	Highly configurable TIS

### iXpress

All the convenience & transportability of a manual wheelchair.  
 All the power & speed of a power chair.



<b>POWER CHAIRS</b>	ADULT									
	Mid Wheel Drive							Rear Wheel Drive		
	Pulse 5	Pulse 6	QM710	QM715 HD	Xperience 2	Xcel 2	QM-710 Stand Up	Salsa R <sup>2</sup>	P-222 SE	Xplore 2
										
<b>Seat Width (mm)</b>	405 - 510	300 - 560	300 - 560	410 - 560	250 - 610	500 - 710	400 - 500	410 - 510	360 - 610	250 - 610
<b>Seat Depth (mm)</b>	405 - 510	300 - 560	300 - 560	410 - 560	250 - 560	430 - 610	400 - 520	410 - 510	360 - 510	250 - 560
<b>Floor to Seat Ht(mm)</b>	400 - 540	400 - 540	400 - 540	440 - 520	400 - 540	420	-	457 - 508	460 - 530	400
<b>Max User Wt (kg)</b>	136	136	136	181	181	249	120	140	161	181
<b>Max Speed (kph)</b>	10	10	10	8	10	8.5	10	10	10	10
<b>Base W x L (mm)</b>	620 x 860	620 x 860	620 x 910	630 x 910	630 x 1100	700 x 1100	630 x 1100	620 x 1140	620 x 800	650 x 1100
<b>Features</b>	VR2: LED single function controls	Multiple control options. Tilt and recline. Power legs	Fully extendable controls. Anti sheer recline and tilt. Suited to outdoors	HD version of QM710	Fully extendable controls. Anti sheer recline and tilt. Low seat to floor height	HD version of the Xperience	Standing wheelchair with superior manoeuvrability. Indoor and outdoor high performance with a multi-fuctional seat. Spidertrac suspension	Fully adjustable seating system. Powered tilt mode. Suited to outdoors	Fully expandable controls. Suits active user. Fixed seat width	Fully extendable controls. Anti sheer recline and tilt. Low seat to floor height

<b>POWER CHAIRS</b>	PAEDIATRIC
	Mid Wheel Drive
	ZM-310
	
<b>Seat Width (mm)</b>	250 - 410
<b>Seat Depth (mm)</b>	250 - 480
<b>Floor to Seat Ht (mm)</b>	390 - 550
<b>Max User Wt (kg)</b>	75
<b>Max Speed (kph)</b>	10
<b>Base W x L (mm)</b>	610 x 890
<b>Features</b>	Power tilt, seat elevate, manual recline, fully expandable controls


  

  

  
