



MANUAL CHAIRS	ADULT								PAEDIATRIC
	Folding		Rigid						
	2GX2	Aero X2	TR3	Aero T	ZRA Series 2	ZR Series 2	Aero Z Series 2	TRA	Twist
Seat Width (mm)	305 - 560	305 - 560	305 - 510	305 - 510	250 - 510	250 - 510	250 - 510	305 - 510	200 - 380
Seat Depth (mm)	305 - 510	305 - 510	305 - 510	305 - 510	250 - 510	250 - 510	250 - 510	305 - 510	200 - 410
Floor to Seat Ht (mm)	340 - 535	340 - 545	380 - 530	380 - 530	405 - 530	405 - 530	405 - 530	380 - 530	320 - 530
Max User Weight (kg)	120	120	120	120	120	120	120	120 / 158 HD	75
Weight (kg)	5.9	7.1	4.2	5.4	4.9	4.2	5.6	5	5.3

Related Products:



SmartDrive



Bodypoint

FreeWheel



the ROHO group
shape fitting technology®

