


Pressure Care Assist



- Can be used anywhere in the bed to assist with pressure care
- Made from medical grade sheepskin
- Helps reduce the risk of pressure area breakouts in the heels, sacrum, shoulders or head

Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_098800044_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity
PTA521600	Pressure Care Assist	156	1