

## Rehband 4143 Manu 3D Basis Elastic Wrist Brace



- High quality 3 tab wrist brace with flexion stay
- Allows free thumb motion
- Six sizes left and right

 Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Wrist Circumference
FBS301385	Rehband 4143 Manu 3D Basis Wrist Brace - Right - X Sml	14.5 - 15.5 cm
FBS301385LL	Rehband 4143 Manu 3D Basis Wrist Brace - Left - Lge	17 - 18 cm
FBS301385LM	Rehband 4143 Manu 3D Basis Wrist Brace - Left - Med	16 - 17 cm
FBS301385LS	Rehband 4143 Manu 3D Basis Wrist Brace - Left - Sml	15.5 - 16 cm
FBS301385LXL	Rehband 4143 Manu 3D Basis Wrist Brace - Left - XL	18 - 19 cm
FBS301385LXS	Rehband 4143 Manu 3D Basis Wrist Brace - Left - X Sml	14.5 - 15.5 cm
FBS301385LXXL	Rehband 4143 Manu 3D Basis Wrist Brace - Left - XXL	19 - 20 cm

<b>Code</b>	<b>Display Name</b>	<b>Wrist Circumference</b>
FBS301385RM	Rehband 4143 Manu 3D Basis Wrist Brace - Right - Med	16 - 17 cm
FBS301385RS	Rehband 4143 Manu 3D Basis Wrist Brace - Right - Sml	15.5 - 16 cm
FBS301385RXL	Rehband 4143 Manu 3D Basis Wrist Brace - Right - XL	18 - 19 cm
FBS301385RXXL	Rehband 4143 Manu 3D Basis Wrist Brace - Right - XXL	19 - 20 cm