

## K Care Shower Commode Weight Bearing Platform



- Assists transportation of the user for showering and can be positioned over the toilet for toileting
- Weight bearing platform footrest and swing out arms assists safe patient transfers in and out of commode
- Ideal for multi user situations such as nursing facilities and hospitals
- Padded or plastic flip up toilet seat available
- Mounted on 4 castors (2 are locking for stability)
- Stainless steel frame for corrosion resistance
- \*Self Propelled model also available
- Note: Pictured with an optional closed padded commode seat (BTC075500)
- **SEATS, PANS & PAN-CARRIERS ARE SOLD AS SEPARATE ITEMS TO THE COMMODE**



### Support @ Home

**Function:** Washing, bathing and showering

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

## Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTC070650	Attendant Propelled Commode - Weight Bearing Platform - 440 mm Wide	150 kg	12 kg	430 mm	440 mm	590 mm
BTC077950	Self Propelled Commode - Weight Bearing Platform - 440 mm Wide	150 kg	17 kg	430 mm	440 mm	590 mm