

## Aspire Watson Chair



The Aspire Watson is an Adjustable Pressure Reducing Chair, designed to deliver exceptional comfort, configurability, and safety across a wide range of care environments. Built with a robust steel frame and industry-leading adjustability, it is an ideal solution for users requiring pressure seating that can be adjusted to suit diverse needs.



- **High-SWL Steel Frame:** Robust steel frame achieves a Safe Working Load (SWL) of 200kg, making it ideal for heavy-duty and long-term care environments.
- **Pressure Injury Prevention:** Supportive memory foam and four-way PU stretch fabric assists with immersion and pressure redistribution during extended use.
- **Optimised Postural Support:** Contoured cushioning combined with adjustable head, seat, and back support promotes proper alignment and sustained comfort.
- **Exceptional Adjustability:** Fully adjustable seat width, depth, and height to accommodate varied body shapes, sizes, and postural needs.
- **Tool-free Adjustment Points:** Ensure quick configuration for multiple users.
- **Integrated Spring-pin Safety System:** Prevents overextension and maintains secure, controlled adjustments.
- **Multi-layered Commercial-grade Foam:** Enhances longevity and sitting comfort.
- **Mobility-ready Rear Castors:** Helps transport the chair between rooms.

### Warranty Information

- 1 Year

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP208040	Aspire Watson Chair - Ash Grey	690 - 790mm	1160 - 1260mm	200kg	Ash Grey	400 - 570mm	500 - 600mm	460 - 570mm	460 - 570mm	800mm