

Oppo 1230 Adjustable Knee Stabilizer



- Two medial and lateral spiral stays provide additional support and stability for the knee
- The additional circle pad around the patella provides padded protection and stabilization to prevent patella displacement
- Open patella design to relieve pressure
- Heat retention and compression reduces potential for strain or injury

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS301524	Oppo - Knee Brace - Coolprene Adj Wrap w/ Stays - O/Patella (UNIVERSAL)	318 - 406 mm