

HipSaver Interim 3/4 Length Overpant



- HipSaver Interim is the perfect solution for hip protection as you implement each patient's individual falls protection programme
- 3/4 length unisex pants with sewn-in hip protection pads
- Immediate hip protection for new admissions, acute care and short stay patients
- As a 'stop gap' component of your Falls Safety Programme while you order regular hipsavers for ongoing hip protection
- HipSaver Interim provides is ideal for:
 - New admissions
 - Acute care and short stay patients
 - Residents with temporary changes in risk factors

- ✔ Contracted DVA product
- ✔ Contract approved

Warranty Information

- 1 Year

Code	Display Name	Waist/Hip Circumference
PTA534107	HipSaver - Interim 3/4 Length Overpant - X Sml	71-81cm
PTA534107S	HipSaver - Interim 3/4 Length Overpant - Sml	82-92cm
PTA534107M	HipSaver - Interim 3/4 Length Overpant - Med	93-102cm
PTA534107L	HipSaver - Interim 3/4 Length Overpant - Lge	103-114cm
PTA534107XL	HipSaver - Interim 3/4 Length Overpant - XL	115-127cm

Code	Display Name	Waist/Hip Circumference
PTA534107XXL	HipSaver - Interim 3/4 Length Overpant - XXL	128-146cm
PTA534108	HipSaver - Interim 3/4 Length Overpant w/ Tailbone Pad - X Sml	71-81cm
PTA534108S	HipSaver - Interim 3/4 Length Overpant w/ Tailbone Pad - Sml	82-92cm
PTA534108M	HipSaver - Interim 3/4 Length Overpant w/ Tailbone Pad - Med	93-102cm
PTA534108L	HipSaver - Interim 3/4 Length Overpant w/ Tailbone Pad - Lge	103-114cm
PTA534108XL	HipSaver - Interim 3/4 Length Overpant w/ Tailbone Pad - XL	115-127cm