

Exercise Bike



- 8 Levels Manual Resistance
- 4 kg Flywheel
- 1 Piece Crank
- 100kg Max User Weight
- Adjustable Stabilisers for uneven surfaces
- Large Cushioned Seat for Comfort
- Vertical Adjustable Seat
- Wheels for Easy Transportation
- LCD Interface
- Hand Pulse Sensors
- Computer Displays: Scan, Time, Speed, Distance, Calories, Pulse, ODO

✓ Contracted DVA product



Compliance

No applicable compliance standard

| Code | Display Name | RRP | Box Quantity | Colour |
|-----------|--------------------------------|-----|--------------|--------|
| FPF290610 | Resistance Level Exercise Bike | 559 | 1 | Black |