

Exercise Bike



- 8 Levels Manual Resistance
- 4 kg Flywheel
- 1 Piece Crank
- 100kg Max User Weight
- Adjustable Stabilisers for uneven surfaces
- Large Cushioned Seat for Comfort
- Vertical Adjustable Seat
- Wheels for Easy Transportation
- LCD Interface
- Hand Pulse Sensors
- Computer Displays: Scan, Time, Speed, Distance, Calories, Pulse, ODO

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Colour
FPF290610	Resistance Level Exercise Bike	Black