

## DermaSaver Forearm Tube Protection



- Soft and loose with ample air flow reaching the skin allowing for greater healing ability whilst offering protection
- Two cushion layers of MicroSpring Textile protect fragile skin from friction and pressure that can cause skin tears, bruising and abrasions
- Sizing: simply measure the circumference around the widest part of the forearm. For borderline measurements, choose the larger size to ensure comfort

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533306	DermaSaver - Forearm Tube - X Sml	13-18cm
PTA537140	DermaSaver - Arm Bow Tube - X Sml	-
PTA533306S	DermaSaver - Forearm Tube - Sml	18-25cm
PTA537100	DermaSaver - Arm Bow Tube - Sml	-
PTA533306M	DermaSaver - Forearm Tube - Med	25-30cm
PTA537110	DermaSaver - Arm Bow Tube - Med	-
PTA533306L	DermaSaver - Forearm Tube - Lge	30-38cm

Code	Display Name	Calf/Arm Circumference
PTA537120	DermaSaver - Arm Bow Tube - Lge	-
PTA533306XL	DermaSaver - Forearm Tube - XL	38-46cm
PTA537130	DermaSaver - Arm Bow Tube - XL	-