

## Orbi-Turn Patient Turner



- The Orbi-turn has been designed for carers to facilitate standing and seat-to-seat transfers
- It allows a patient with moderate upper and lower limb strength to participate in functional transfers whilst ensuring good working postures for carers and therapists
- Safe Working Load 190 kg

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Changing body position

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123106154\_0105\_1\_2

**Minor Transfer –**

**Composite:**05\_430\_0103\_1\_2

Click [here](#) for more information

Code	Display Name	Material
TSE672935	Orbi-Turn Standing Patient Turner	Plastic