

Solaris ReadyWrap - Knee



- The Solaris Lymphoedema ReadyWrap Knee Wrap features three strap design that provides a comfortable fit reducing pinching and high-pressure points
- Provide higher working and lower resting compression levels with low-stretch materials
- Provides support with localized patella swelling, or acts as a bridge connecting the calf and thigh units to provide full leg coverage
- Designed to provide even and effective compression
- Hook & loop strapping for a customized fit. Easily adjustable to allow for daily girth fluctuations

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
FBS303629	Solaris ReadyWrap - Knee - Large - Beige	Beige
FBS303630	Solaris ReadyWrap - Knee - Large - Black	Black
FBS303631	Solaris ReadyWrap - Knee - Medium - Beige	Beige
FBS303632	Solaris ReadyWrap - Knee - Medium - Black	Black
FBS303633	Solaris ReadyWrap - Knee - Small - Beige	Beige
FBS303634	Solaris ReadyWrap - Knee - Small - Black	Black
FBS303635	Solaris ReadyWrap - Knee - Xlarge - Beige	Beige

Code	Display Name	Colour
FBS303637	Solaris ReadyWrap - Knee - XXlarge - Beige	Beige
FBS303638	Solaris ReadyWrap - Knee - XXlarge - Black	Black