

## My Comfort Chair



**Air comfort**

- Tilt / Recline / Elevating Leg-rest
- Visco elastic memory foam throughout minimises pressure without compromising postural support
- Careflex three-way stretch fabric minimises sheer (friction) whilst promoting air-flow
- The Anti-Microbial, Anti-Bacterial properties of the fabric help with infection control and skin health
- Tilt-in Space adjustable seat positioning with infinite locking
- Adjustable wings and arms for easy transfer
- Polyurethane, Anti-Microbial, Anti- Bacterial, Three-way stretch fabric
- Suits users: 162 cm to 188 cm
- Available in Small - Suits users up to 142 cm
- Adjustable backrest for full recline
- Adjustable legrest with infinite raise and lower positions helps position the user to help minimise fluid retention
- Four low-friction locking castors -Three braking, One directional
- Soft-touch brake and recline levers for easier operation

✓ Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

## Warranty Information

- 3 Years Frame
- 1 Year Fabric

Code	Display Name	Safe Working Load	Material	Seat Width	Backrest Height
CHP198390	My Comfort Chair	180 kg	Steel Frame, Careflex	510 mm	710 mm

Code	Display Name	Safe Working Load	Material	Seat Width	Backrest Height
------	--------------	-------------------	----------	------------	-----------------