

Mangar Camel - Lift Cushion



The Camel Lifting Cushion will both sit up and lift a fallen person. Individuals with cognitive impairments find comfort in the Camel's gentle lifting mechanism, while the added backrest support caters perfectly to bariatric patients. It is not just about assistance – it is about elevating the quality of life for everyone involved.

- The lifting cushion is inflated using an Airflo Compressor with a simple hand control. No extensive training is required – it is designed for hassle-free operation, whether you are using it independently or with assistance
- Perfect for seamless transitions between indoor and outdoor settings, it is all about convenience at your fingertips.
- Easy to clean using a proprietary, non-abrasive mild liquid cleaner, simply deflate, fold, and stow away in the provided carry bag – ideal for infection control management.
- With a lifting capacity of up to 320kg, it ensures the safe lifting of even the most fragile individuals from the floor. Its stable and supportive design is a boon for caregivers dealing with individuals prone to falls, offering unmatched peace of mind.
- Experience a lift in both safety and dignity. By reducing the risk of musculoskeletal injuries for caregivers, it's a game-changer in care routines.

✓ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name
------	--------------

LSS403810	Mangar Camel - Lift Cushion - 320kg
-----------	-------------------------------------