

Leckey Squiggles Standing System



- Squiggles Standing System supports the child in a comfortable, secure, upright, prone or supine standing posture
- Chassis can clip on and off in seconds for folding and storage
- Adjustable to accommodate children with differing needs
- Accessories available include tray, knee blocks, supports, headrests which all attach without tools
- User Standing Height: 750 - 1100 mm
- Chest Pad Height: 550 - 820 mm



ndis

Item Code:

05_053603010_0103_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

LECKEY[®]

Code	Display Name	Safe Working Load
PED479379	Leckey Squiggles Standing System, Ages 1-5 years	22 kg