

## DermaSaver Arm Tube with Knuckle Protector



- Soft and loose with ample airflow reaching the skin allowing for greater healing ability whilst offering protection
- Two cushion layers of MicroSpring Textile protect fragile skin from friction and pressure that can cause skin tears, bruising and abrasions, plus added knuckle protection
- Sizing: simply measure the circumference around the widest part of the forearm. For borderline measurements, choose the larger size to ensure comfort

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533311	DermaSaver - Arm Tube with Knuckle Protector - X Sml	13-18cm
PTA533311S	DermaSaver - Arm Tube with Knuckle Protector - Sml	18-25cm
PTA533311M	DermaSaver - Arm Tube with Knuckle Protector - Med	25-30cm
PTA533311L	DermaSaver - Arm Tube with Knuckle Protector - Lge	30-38cm
PTA533311XL	DermaSaver - Arm Tube with Knuckle Protector - XL	38-46cm