

Alerta - Wireless - Nurse Call Button System



- Wireless nurse call system is a kit that includes call button and alarm monitor
- Provides the user with reassurance that help is at hand. With a large easy to press button it activates the alarm when a user presses the button. This alerts the caregiver to ascertain the user's need
- Connects to a nurse call system through the Alerta Alarm Monitor for centralised activation and comes with a lanyard, and a holder that can be wall mounted, so it can always be used remotely within the homecare environment
- Simple plug and play ensures carer can quickly set up and operate with ease
- High and low volume alert

✔ Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year



Code

Display Name

BEA007631 Alerta - Wireless - Nurse Call Button System (Includes button & receiver)