

MoliCare Premium Elastic



All-in-one incontinence slip with 2 elasticated side panels: suits and adapts to all body shapes and movements and provides an ergonomic use - approximately 20% faster application as well as easier and more ergonomic handling for caregivers.

- Comfortable, adjusted fit and ergonomic handling due to the 2 elasticated side panels and re-closable tapes
- Skin-friendly thanks to soft and breathable material
- The antibacterial finish and pH buffering help to maintain a skin-neutral pH value of 5.5
- All-round protection is provided by fast absorption, anti-leakage system and quick dry system
- Hypoallergenic and dermatologically tested
- Wetness indicator on the back sheet shows when it is time to change



Compliance

No applicable compliance standard

Code	Display Name	Capacity
CON254300	MoliCare Premium Elastic 8D - Large (Pkt 24)	8 Drops