

## NuboMed M3100-T Medication Management Cabinet



- **User Interaction:** Features a 15.6" touch screen monitor with a Capacitive Fingerprint Reader, optional facial recognition (with mask recognition capability), Staff RFID card reader, Built-in monitor speaker, and a Barcode Scanner (1D/2D barcodes)
- **Cabinet Material:** Constructed from high-strength metal steel with medical white powder coating that has disinfectant-resistant, anti-aging and anti-discoloration properties
- **Cabinet Control Lock:** Intelligent Lock Control; Authorised user access only, Hierarchical Access Control
- **Drawer Layers:** Configurable upon order
- **Proximity Scanner:**
  - Recognition distance: <50cm
  - Recognition speed: <2 seconds
- **Guiding Light:** Each bin has a RGB indicator light as a guide to direct operation
- **Smart Colour Bin Display:** The mini colour screen displays the name, specification, and instructions if available
- **Temperature and Humidity Display:** Features a digitalised LED panel display on the cabinet
- **Temperature and Humidity Alarm:** Smart detection of temperature and humidity levels with an auto-alarm
- **Ideal For:** Outpatient Pharmacy, Emergency Room, Pharmacy, Central Pharmacy and Clinical Ward

### Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
WAT908562	NuboMed M3100-T Medication Management Cabinet	528 mm	529 mm	771 mm	120 kg