

Stealth Spectrum Gel



- Lateral and pre-ischial shelf contoured, high-density moulded foam base
- CoolCore® technology cover standard
- Cushion height is 3"
- Safe working load 226 kg

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTC623100	Stealth Spectrum Gel - 14in x 18in	18 in	14 in	3 in	226 kg	Medium Risk
PTC599801	Stealth Spectrum Gel - 16 x 16in	16 in	16 in	-	-	High Risk
PTC599901	Stealth Spectrum Gel - 16 x 18in	16 in	18 in	-	-	High Risk
PTC623130	Stealth Spectrum Gel - 16in x 19in	19 in	16 in	3 in	226 kg	Medium Risk
PTC600101	Stealth Spectrum Gel - 18 x 16in	18 in	16 in	-	-	High Risk
PTC600201	Stealth Spectrum Gel - 18 x 18in	18 in	18 in	-	-	High Risk
PTC623150	Stealth Spectrum Gel - 17in x 16in	16 in	17 in	3 in	226 kg	Medium Risk

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTC600251	Stealth Spectrum Gel - 18 x 20in	18 in	20 in	-	-	High Risk
PTC623160	Stealth Spectrum Gel - 17in x 17in	17 in	17 in	3 in	226 kg	Medium Risk
PTC600301	Stealth Spectrum Gel - 20 x 22in	20 in	22 in	-	-	High Risk
PTC623170	Stealth Spectrum Gel - 17in x 18in	18 in	17 in	3 in	226 kg	Medium Risk
PTC623180	Stealth Spectrum Gel - 17in x 19in	19 in	17 in	3 in	226 kg	Medium Risk
PTC623190	Stealth Spectrum Gel - 17in x 20in	20 in	17 in	3 in	226 kg	Medium Risk
PTC623210	Stealth Spectrum Gel - 18in x 17in	17 in	18 in	3 in	226 kg	Medium Risk
PTC623240	Stealth Spectrum Gel - 18in x 22in	22 in	18 in	3 in	226 kg	Medium Risk
PTC623250	Stealth Spectrum Gel - 18in x 24in	24 in	18 in	3 in	226 kg	Medium Risk
PTC623260	Stealth Spectrum Gel - 19in x 17in	17 in	19 in	3 in	226 kg	Medium Risk
PTC623270	Stealth Spectrum Gel - 19in x 18in	18 in	19 in	3 in	226 kg	Medium Risk
PTC623280	Stealth Spectrum Gel - 19in x 19in	19 in	19 in	3 in	226 kg	Medium Risk
PTC623290	Stealth Spectrum Gel - 19in x 20in	20 in	19 in	3 in	226 kg	Medium Risk
PTC623300	Stealth Spectrum Gel - 20in x 16in	16 in	20 in	3 in	226 kg	Medium Risk
PTC623310	Stealth Spectrum Gel - 20in x 18in	18 in	20 in	3 in	226 kg	Medium Risk
PTC623320	Stealth Spectrum Gel - 20in x 19in	19 in	20 in	3 in	226 kg	Medium Risk
PTC623330	Stealth Spectrum Gel - 20in x 20in	20 in	20 in	3 in	226 kg	Medium Risk

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTC623350	Stealth Spectrum Gel - 20in x 24in	24 in	20 in	3 in	226 kg	Medium Risk
PTC623360	Stealth Spectrum Gel - 22in x 16in	16 in	22 in	3 in	226 kg	Medium Risk
PTC623370	Stealth Spectrum Gel - 22in x 18in	18 in	22 in	3 in	226 kg	Medium Risk
PTC623380	Stealth Spectrum Gel - 22in x 19in	19 in	22 in	3 in	226 kg	Medium Risk
PTC623400	Stealth Spectrum Gel - 22in x 22in	22 in	22 in	3 in	226 kg	Medium Risk
PTC623410	Stealth Spectrum Gel - 22in x 24in	24 in	22 in	3 in	226 kg	Medium Risk
PTC623420	Stealth Spectrum Gel - 24in x 16in	16 in	24 in	3 in	226 kg	Medium Risk
PTC623430	Stealth Spectrum Gel - 24in x 18in	18 in	24 in	3 in	226 kg	Medium Risk
PTC623440	Stealth Spectrum Gel - 24in x 20in	20 in	24 in	3 in	226 kg	Medium Risk
PTC623450	Stealth Spectrum Gel - 24in x 22in	22 in	24 in	3 in	226 kg	Medium Risk
PTC623460	Stealth Spectrum Gel - 24in x 24in	24 in	24 in	3 in	226 kg	Medium Risk
PTC623255	Stealth Spectrum Gel - 24 x 19in	24 in	19 in	3 in	226 kg	Medium Risk
PTC623115	Stealth Spectrum Gel - 16 x 17in	16 in	17 in	3 in	226 kg	Medium Risk
PTC623215	Stealth Spectrum Gel - 18 x 19in	18 in	19 in	3 in	226 kg	Medium Risk
PTC623455	Stealth Spectrum Gel - 24 x 21in	24 in	21 in	3 in	226 kg	Medium Risk