

TENA Pads Comfort



- Provides an excellent level of security, comfort and discretion
- Best suited to ambulant people or those who are able to toilet themselves
- To be used in conjunction with TENA Fix Stretch Pants or TENA Fix Cotton Special as part of a two-piece system

Warranty Information

- 1 Year



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Capacity
CON243000	TENA Pads Comfort Plus (Pkt 20)	780 ml
CON243100	TENA Pads Comfort Extra (Pkt 20)	1080 ml
CON243200	TENA Pads Comfort Super (Pkt 20)	1320 ml
CON243300	TENA Pads Comfort Maxi (Pkt 20)	1560 ml